



JIMY

DENTAL INSTRUMENTS




PRODUCTS CATALOGUE

Tooth Extracting Forceps|(eng)




Incisors and Cuspids 1
JU-51-001
 Tooth Ext Forceps
 Fig.1,



Incisors and Bicuspids 2
JU-51-002
 Tooth Ext Forceps
 Fig.2,




Lower Incisors and Canines 4
JU-51-004
 Tooth Ext Forceps
 Fig.4,



Bicuspids 7
JU-51-007
 Tooth Ext Forceps
 Fig.7,




Lower Bicuspids, Universal 13
JU-51-013
 Tooth Ext Forceps
 Fig.13,




Molars, Right 17
JU-51-017
 Tooth Ext Forceps
 Fig.17,

Tooth Extracting Forceps|(eng)



Molars, Left 18
JU-51-018
 Tooth Ext Forceps
 Fig.18,



Molars 18A
JU-51-018-01
 Tooth Ext Forceps
 Fig.18A,



Lower Molars, Universal 22
JU-51-022
 Tooth Ext Forceps
 Fig.22,



Lower Molars, Right 23
JU-51-023
 Tooth Ext Forceps
 Fig.23,



Lower Molars, Left 24
JU-51-024
 Tooth Ext Forceps
 Fig.24,



Upper roots and Incisors 29
JU-51-029
 Tooth Ext Forceps
 Fig.29,

Tooth Extracting Forceps|(eng)


Lower Roots 33
JU-51-033
 Tooth Ext Forceps
 Fig.33,



Lower Roots 01
JU-51-033-01
 Tooth Ext Forceps
 Fig.33 A,




Lower Roots 03
JU-51-033-03
 Tooth Ext Forceps
 Fig.33 C,



Lower Roots 06
JU-51-033-06
 Tooth Ext Forceps
 Fig.33 M,



Incisors and Bicuspid 34N
JU-51-034
 Tooth Ext Forceps
 Fig.34,



Incisors and Bicuspid 35N
JU-51-035
 Tooth Ext Forceps
 Fig.35,

Tooth Extracting Forceps|(eng)



JU-51-036
Tooth Ext Forceps
Fig.36,



Lower roots 46
JU-51-046
Tooth Ext Forceps
Fig.46,



Lower Roots 46
JU-51-046-01
Tooth Ext Forceps
Fig.46 L,



Lower Roots 46
JU-51-046-02
Tooth Ext Forceps
Fig.46 LX,



Very Fine Roots 49
JU-51-049
Tooth Ext Forceps
Fig.49,



Roots 51
JU-51-051
Tooth Ext Forceps
Fig.51,

Tooth Extracting Forceps|(eng)



Roots 51A
JU-51-051-01
 Tooth Ext Forceps
 Fig.51 A,



Roots 51 C
JU-51-051-03
 Tooth Ext Forceps
 Fig.51 C,



JU-51-057
 Tooth Ext Forceps
 Fig.57,



JU-51-057-1
 Tooth Ext Forceps
 Fig.57 A,




JU-51-057-02
 Tooth Ext Forceps
 Fig.57 B,



Third Molars 67
JU-51-067
 Tooth Ext Forceps
 Fig.67,

Tooth Extracting Forceps|(eng)



Third Molars 67 A
JU-51-067-01
 Tooth Ext Forceps
 Fig.67 A,



Third Molars Deep Griping 67 N
JU-51-067-02
 Tooth Ext Forceps
 Fig.67 N,



Lower Molars 73
JU-51-073
 Tooth Ext Forceps
 Fig.73,



Lower Roots 74
JU-51-074
 Tooth Ext Forceps
 Fig.74,



Incisors and Roots 74N
JU-51-074-06
 Tooth Ext Forceps
 Fig.74 N,



Incisors and Roots ,Slender Jaw74XN
JU-51-074-07
 Tooth Ext Forceps
 Fig.74 XN,

Tooth Extracting Forceps|(eng)



Incisors and Bicuspid
JU-51-075
 Tooth Ext Forceps
 Fig.75,



Upper Roots 76
JU-51-076
 Tooth Ext Forceps
 Fig.76,




Upper Roots 76N
JU-51-076-07
 Tooth Ext Forceps
 Fig.76 N,



Third Molars 79
JU-51-079
 Tooth Ext Forceps
 Fig.79,



Third Molars 79A
JU-51-079-01
 Tooth Ext Forceps
 Fig.79 A,



Third Molars ,Deep Gripping 79N
JU-51-079-02
 Tooth Ext Forceps
 Fig.79 N,

Tooth Extracting Forceps|(eng)



Molars with carious or broken caps 86A

JU-51-086-01

Tooth Ext Forceps
Fig.86 A,



Molars with carious or broken caps 86B

JU-51-086-02

Tooth Ext Forceps
Fig.86 B,



Molars with carious or broken caps 86C

JU-51-086-03

Tooth Ext Forceps
Fig.86 C,



JU-51-097

Tooth Ext Forceps
Fig.97,



Very Fine Roots 133

JU-51-133

Tooth Ext Forceps
Fig.133,



Roots Fragments, Deep Gripping 146

JU-51-146

Tooth Ext Forceps
Fig.146,

Tooth Extracting Forceps|(eng)



Tooth Ext Forceps
JU-51-151
 Tooth Ext Forceps
 Fig.151,



Tooth Ext Forceps
JU-51-153
 Fig.153,



Roots, with Parallel beak 233
JU-51-233
 Tooth Ext Forceps
 Fig.233,



Incisors and Bicuspid 1
JU-51-301
 Tooth Ext Forceps
 Fig.1,




Incisors and Bicuspid 2
JU-51-302
 Tooth Ext Forceps
 Fig.2,



Bicuspid 7
JU-51-307
 Tooth Ext Forceps
 Fig.7,

Tooth Extracting Forceps|(eng)


Bicuspid 13
JU-51-313
 Tooth Ext Forceps
 Fig.13,



Molars 17
JU-51-317
 Tooth Ext Forceps
 Fig.17,



Molars 18
JU-51-318
 Tooth Ext Forceps
 Fig.18,



Molars 22
JU-51-322
 Tooth Ext Forceps
 Fig.22,



Roots 33
JU-51-330
 Tooth Ext Forceps
 Fig.30,



Roots 33
JU-51-333
 Tooth Ext Forceps
 Fig.33,

Tooth Extracting Forceps|(eng)


Roots 51
JU-51-351
 Tooth Ext Forceps
 Fig.51,



Roots, Very Fine Beak 51A
JU-51-351-01
 Tooth Ext Forceps
 Fig.51 A,



Third Molars 67 A
JU-51-367-01
 Tooth Ext Forceps
 Fig.67 A,



Molars 73
JU-51-373
 Tooth Ext Forceps
 Fig.73,




Incisors and Roots 74
JU-51-374
 Tooth Ext Forceps
 Fig.74,



Incisors and roots , with fine beak 74N
JU-51-374-01
 Tooth Ext Forceps
 Fig.74 N,

Tooth Extracting Forceps|(eng)



Third Molars 79

JU-51-379

Tooth Ext Forceps
Fig.79,



Incisors and cuspids 1

JU-51-401

Tooth Ext Forceps
Fig.1 Diam,



Bicuspid 7

JU-51-407

Tooth Ext Forceps
Fig.7 Diam,



Incisors and cuspids 13

JU-51-413

Tooth Ext Forceps
Fig.13 Diam,



Molars , Right 17

JU-51-417

Tooth Ext Forceps
Fig.17 Diam,



Molars , left 18

JU-51-418

Tooth Ext Forceps
Fig.18 Diam,

Tooth Extracting Forceps|(eng)




Molars 22
JU-51-422
 Tooth Ext Forceps
 Fig.22 Diam,



Roots 33
JU-51-433
 Tooth Ext Forceps
 Fig.33 Diam,



Very Fine Roots 46
JU-51-446
 Tooth Ext Forceps
 Fig.46 Diam,



Very Fine Roots 49
JU-51-449
 Tooth Ext Forceps
 Fig.49 Diam,



Roots 51
JU-51-451
 Tooth Ext Forceps
 Fig.51 Diam,



Roots , With Fine Beak 51A
JU-51-451-01
 Tooth Ext Forceps
 Fig.51 A Diam,

Tooth Extracting Forceps|(eng)




Martin Ergo 67A
JU-51-467-01
 Tooth Ext Forceps
 Fig.67 A Diam,



Martin Ergo 74N
JU-51-474
 Tooth Ext Forceps
 Fig.74 N Diam,



Martin Ergo 79
JU-51-479
 Tooth Ext Forceps
 Fig.79 Diam,



Martin Ergo ,Very Fine Roots 97
JU-51-497
 Tooth Ext Forceps
 Fig.97 Diam,



Bicuspid 13S
JU-51-613-09
 Tooth Ext Forceps
 Fig.13 S,



Molars 22s
JU-51-622-09
 Tooth Ext Forceps
 Fig.22 S,

Tooth Extracting Forceps|(eng)



Incisors and Roots 30S
JU-51-630-09
 Tooth Ext Forceps
 Fig.30 S,



Roots 33S
JU-51-633-09
 Tooth Ext Forceps
 Fig.33 S,



Incisors and bicuspid 37
JU-51-637
 Tooth Ext Forceps
 Fig.37,



Incisors and Cuspids 38
JU-51-638
 Tooth Ext Forceps
 Fig.38,



Molars 39
JU-51-639
 Tooth Ext Forceps
 Fig.39,



Molars 39A
JU-51-639-01
 Tooth Ext Forceps
 Fig.39 A,

Tooth Extracting Forceps|(eng)



Molars , Left 39L
JU-51-639-05
 Tooth Ext Forceps
 Fig.39 L,



Molars , Right 39R
JU-51-639-08
 Tooth Ext Forceps
 Fig.39 R,



Molars 40
JU-51-640
 Tooth Ext Forceps
 Fig.40,



Roots 51S
JU-51-651-09
 Tooth Ext Forceps
 Fig.51 S,



Roots 76S
JU-51-676-09
 Tooth Ext Forceps
 Fig.76 S,



Molars 3
JU-51-703
 Tooth Ext Forceps
 Fig.3,

Tooth Extracting Forceps|(eng)



Incisors 5
JU-51-705
 Tooth Ext Forceps
 Fig.5,



Molars 6
JU-51-706
 Tooth Ext Forceps
 Fig.6,



Roots 7
JU-51-707
 Tooth Ext Forceps
 Fig.7,



Incisors and cuspids 137
JU-51-737
 Tooth Ext Forceps
 Fig.137,



Bicuspid 139
JU-51-739
 Tooth Ext Forceps
 Fig.139,



Klein , Roots 51S
JU-51-751-09
 Tooth Ext Forceps
 Fig.51 S,

Tooth Extracting Forceps|(eng)



Molars 17SK
JU-51-817-09
Tooth Ext Forceps
Fig.17 SK,



Molars 23SK
JU-51-823-09
Tooth Ext Forceps
Fig.23 SK,



Teeth and Roots ,Universal 150SK
JU-51-850-09
Tooth Ext Forceps
Fig.150 SK,



Teeth and Roots, Universal 151Sk
JU-51-851-09
Tooth Ext Forceps
Fig.151 SK,

Tooth Extracting Forceps|(amr)


Incisors and Bicuspid 1
VI-52-001
 Tooth Ext Forceps Amr
 Fig.1,



Henahan Incisors 1A
VI-52-001-01
 Tooth Ext Forceps Amr
 Fig.1 A,



Molars and third molar 10S
JU-52-010-09
 Tooth Ext Forceps Amr
 Fig.10 S,



Molars 15
JU-52-015
 Tooth Ext Forceps Amr
 Fig.15,



Molars 16
JU-52-016
 Tooth Ext Forceps Amr
 Fig.16,



Molars , Slender Jaw
JU-52-016-09
 Tooth Ext Forceps Amr
 Fig.16 S,

Tooth Extracting Forceps|(amr)


Molars 17
JU-52-017
 Tooth Ext Forceps Amr
 Fig.17,



Harris Molars , left 18L
JU-52-018-05
 Tooth Ext Forceps Amr
 Fig.18 L,



Harris Molars , Right 18R
JU-52-018-08
 Tooth Ext Forceps Amr
 Fig.18 R,



Molars
JU-52-023
 Tooth Ext Forceps Amr
 Fig.23,




Parmlly Molars and bicuspid 32
JU-52-032
 Tooth Ext Forceps Amr
 Fig.32,




Parmlly Molars and Bicuspid , Roots 32A
JU-52-032-01
 Tooth Ext Forceps Amr
 Fig. 32 A,

Tooth Extracting Forceps|(amr)



Molars , Left 53L
JU-52-053-05
 Tooth Ext Forceps Amr
 Fig.53 L,



Molars , Right 53R
JU-52-053-08
 Tooth Ext Forceps Amr
 Fig.53 R,



Teeth and Roots, Universal also for Children 62
JU-52-062
 Tooth Ext Forceps Amr
 Fig.62,



Incisors and Roots 65
JU-52-065
 Tooth Ext Forceps Amr
 Fig.65,



Tomes Roots Fragments and Small Roots 69
JU-52-069
 Tooth Ext Forceps Amr
 Fig.69,



Nevius Molars , Left 88L
JU-52-088-05
 Tooth Ext Forceps Amr
 Fig.88 L,

Tooth Extracting Forceps|(amr)


Nevius Molars , Right 88R
JU-52-088-08
 Tooth Ext Forceps Amr
 Fig.88 R,



Cryer Bicuspid, Incisors and Roots 150
JU-52-150
 Tooth Ext Forceps Amr
 Fig.150,



Cryer Bicuspid, Incisors and Roots 150A
JU-52-150-01
 Tooth Ext Forceps Amr
 Fig.150 A,



Bicuspid and Incisors , Roots 150X
JU-52-150-02
 Tooth Ext Forceps Amr
 Fig.150 X,



Bicuspid and Incisors, Roots 150AS
JU-52-150-08
 Tooth Ext Forceps Amr
 Fig.150 AS,



Cryer Bicuspid and Incisors , Roots 150S
JU-52-150-09
 Tooth Ext Forceps Amr
 Fig.150 S,

Tooth Extracting Forceps|(amr)



Cryer Bicuspid , Incisors and Roots 151

JU-52-151

Tooth Ext Forceps Amr
Fig.151,



Cryer Bicuspid , Incisors and Roots 151A

JU-52-151-01

Tooth Ext Forceps Amr
Fig.151 A,



Cryer Bicuspid , Incisors and Roots 151X

JU-52-151-02

Tooth Ext Forceps Amr
Fig.151 X,



Bicuspid and Incisors , Roots 151AS

JU-52-151-08

Tooth Ext Forceps Amr
Fig.151 AS,



Cryer Teeth and Roots , Universal 151S

JU-52-151-09

Tooth Ext Forceps Amr
Fig.151 S,



Third Molars 222

JU-52-222

Tooth Ext Forceps Amr
Fig.222,

Tooth Extracting Forceps|(amr)



Root Fragments and Universal 501

JU-52-501

Tooth Ext Forceps Amr
Fig.501,



Witzel Root Fragments and Universal 502

JU-52-502

Tooth Ext Forceps Amr
Fig.502,

Root Elevators



Seldin 1L
JU-53-001-05
 Root Elevators
 Fig.IL,




Seldin 1R
JU-53-001-08
 Root Elevator
 Fig.IR,



Seldin 2
JU-53-002-00
 Root Elevators
 Fig.2,



Seldin
JU-53-003-00
 Root Elevators
 Fig.3,



Cryer
JU-53-003-01
 Root Elevators
 Fig.3 A,



Seldin 4L
JU-53-004-05
 Root Elevators
 Fig. 4 L,

Root Elevators



Seldin 4R
JU-53-004-08
 Root Elevators
 Fig. 4 R,



Crane Pick 6
JU-53-006-00
 Root Elevators
 Fig. 6,



Miller 1
JU-53-008-00
 Root Elevators
 Fig. 8,



Miller 2
JU-53-009-00
 Root Elevators
 Fig.9,



Miller-Apexo
JU-53-010-00
 Root Elevators
 Fig.10,



Miller-Apexo
JU-53-011-00
 Root Elevators
 Fig.11,

Root Elevators


Cryer
JU-53-012-00
 Root Elevators
 Fig.12,




Cryer
JU-53-027-00
 Root Elevators
 Fig.27,



Cryer 28
JU-53-028-00
 Root Elevators
 Fig.28,



Seldin 34
JU-53-034-00
 Root Elevators
 Fig.34,



Seldin 34S
JU-53-034-09
 Root Elevators
 Fig.34 S,




Cryer 27S
JU-53-039-00
 Root Elevators
 Fig. 39,

Root Elevators



Cryer 28S
JU-53-040-00
 Root Elevators
 Fig. 40,



Flohr
JU-53-046-00
 Root Elevators
 Fig. 46,



Barry 50
JU-53-050-00
 Root Elevators
 Fig. 50,



Barry 51
JU-53-051-00
 Root Elevators
 Fig. 51,



Lindo Levian 1
JU-53-200-01
 Root Elevators
 Fig. 1,




Lindo Levian 3
JU-53-200-03
 Fig. 3,


Root Elevators




Howard-Apexo # 1
JU-53-204-00
 Root Elevators
 Fig. 1,



JU-53-204-01
 Root Elevators
 Fig. 1,




Root Elevators
JU-53-205-00
 Fig. 1,



Root Elevators
JU-53-205-01
 Fig. 1,




Bein 1
JU-53-220-01
 Root Elevator
 fig. 2,




Bein 2
JU-53-220-02
 Root Elevator
 Fig. 3,


Root Elevators




Bein 3
JU-53-220-03
 Root Elevator
 Fig. 3,




Bein 4
JU-53-220-04
 Fig. 4,




Bein 1A
JU-53-221-01
 Fig. 1 A,



Bein 2A
JU-53-221-02
 Fig. 2 A,




Bein 3A
JU-53-221-03
 Fig. 3 A,




Bein 4A
JU-53-221-04
 Fig. 4 A,

Root Elevators


Apical 1
JU-53-222-01
Fig. 1,




Apical 2
JU-53-222-02
Fig. 2,




Apical 3
JU-53-222-03
Fig. 3,



Apical 4
JU-53-222-04
Fig. 4,




Apical 1A
JU-53-223-01
Fig. 1 A,



Apical 2A
JU-53-223-02
Fig. 2 A,

Root Elevators



Apical 3A
JU-53-223-03
 Fig. 3 A,



Apical 4A
JU-53-223-04
 Fig. 4 A,




Apical 301
JU-53-301-00
 Root Elevators
 Fig. 301,




Apical 302
JU-53-302-00
 Root Elevators
 Fig. 302,



Apical-Flohr
JU-53-303-00
 Root Elevators
 Fig. 303,



Apical 304
JU-53-304-00
 Root Elevators
 Fig. 304,

Root Elevators


Apical 46R
JU-53-311-00
 Root Elevators
 Fig. 46 R,



Coleman
JU-53-405-01
 Root Elevators
 Fig. 1,




Coleman 2L
JU-53-405-02
 Root Elevators
 Fig. 2 L,




Coleman 3R
JU-53-405-03
 Root Elevators
 Fig. 3 R,



Coupland Bone chisels 1
JU-53-415-01
 Root Elevators
 Fig. 1,



Coupland Bone Chisels 2
JU-53-415-02
 Root Elevators
 Fig. 2,

Root Elevators


Coupland Bone chisels 3
JU-53-415-03
 Root Elevators
 Fig. 3,



Friedman 1
JU-53-505-01
 Root Elevators
 Fig. 1,




Friedman 2L
JU-53-505-02
 Root Elevators
 Fig. 2 L,




Friedman 2R
JU-53-505-03
 Root Elevators
 Fig. 3 R,




Hylin 1
JU-53-508-01
 Root Elevators
 Fig. 1,



Hylin 2
JU-53-508-02
 Root Elevators
 Fig. 2,

Root Elevators



Hylin 3L
JU-53-508-03
 Root Elevators
 Fig. 3 L,




Hylin 4R
JU-53-508-04
 Root Elevators
 Fig. 4 R,



Bein 0
JU-53-511-00
 Root Elevators
 Fig. 0,



Bein 1
JU-53-511-01
 Root Elevators
 Fig. 1,




Bein 2
JU-53-511-02
 Root Elevators
 Fig. 2,




Bein
JU-53-512-00
 Root Elevators
 Fig. 0,

Root Elevators




Bein 1
JU-53-512-01
 Root Elevators
 Fig. 1,




Bein 2
JU-53-512-02
 Root Elevators
 Fig. 2,



Bein 0
JU-53-513-00
 Root Elevators
 Fig. 0,



Bein 1
JU-53-513-01
 Root Elevators
 Fig. 1,




Bein 2
JU-53-513-02
 Root Elevators
 Fig. 2,




Bein 0A
JU-53-514-00
 Root Elevators
 Fig. 0 A,

Root Elevators




Bein 1A
JU-53-514-01
 Root Elevators
 Fig. 1 A,




Bein 2A
JU-53-514-02
 Root Elevators
 Fig. 2 A,



Bein 0A
JU-53-515-00
 Root Elevators
 Fig. 0 A,



Bein 1A
JU-53-515-01
 Root Elevators
 Fig. 1 A,



Bein 2A
JU-53-515-02
 Root Elevators
 Fig. 2 A,



W.F. Barry 320L
JU-53-533-20
 Root Elevators
 Fig. 320 L,

Root Elevators


W.F. Barry 321R
JU-53-533-21
 Root Elevators
 Fig. 321 R,



Barry 322
JU-53-533-22
 Root Elevators
 Fig. 322,



Barry 323
JU-53-533-23
 Root Elevators
 Fig. 323,



W.F. Barry 370R
JU-53-533-70
 Root Elevators
 Fig. 370 R,



W.F. Barry 371R
JU-53-533-71
 Root Elevators
 Fig. 371 R,



Winter 11L
JU-53-551-05
 Root Elevators
 Fig. 11 L,

Root Elevators



Winter 11R
JU-53-551-08
 Root Elevators
 Fig. 11 R,



Winter 13L
JU-53-553-05
 Root Elevators
 Fig. 13 L,




Winter 13R
JU-53-553-08
 Root Elevators
 Fig. 13 R,



Winter 14L
JU-53-554-05
 Root Elevators
 Fig. 14 L,



Winter 14R
JU-53-554-08
 Root Elevators
 Fig. 14 R,



Lecluse 1
JU-53-561-01
 Root Elevators
 Fig. 1,

Root Elevators


Pott 1
JU-53-567-01
 Root Elevators
 Fig. 1,



Pott 2
JU-53-567-02
 Root Elevators
 Fig. 2,



Krallenheber 1L
JU-53-571-01
 Root Elevators
 Fig. 1 L,



Krallenheber 2R
JU-53-571-02
 Root Elevators
 Fig. 2 R,



Warwick-James 1
JU-53-576-01
 Root Elevators
 Fig. 1,




Warwick-James 2L
JU-53-576-02
 Root Elevators
 Fig. 2 L,

Root Elevators



Warwick-James 3R
JU-53-576-03
 Root Elevators
 Fig. 3 R,



Chompret 1
JU-53-577-01
 Root Elevators
 Fig. 1,



Chompret 2
JU-53-577-02
 Root Elevators
 Fig. 2,



Chompret 3
JU-53-577-03
 Root Elevators
 Fig. 3,




Heidebrink 1
JU-53-582-01
 Root Elevators
 Fig. 1,



Heidebrink 2
JU-53-582-02
 Root Elevators
 Fig. 2,


Root Elevators




Heidebrink 3
JU-53-582-03
Root Elevators
Fig. 3,



West 1
JU-53-584-01
Root Elevators
Fig. 1,



West 2
JU-53-584-02
Root Elevators
Fig. 2,



West 3
JU-53-584-03
Root Elevators
Fig. 3,

Bone Rongeurs


Bone Rongeur
JU-54-101
 Bone Holding Forceps
 185 mm, 7?",




Bone Rongeur
JU-54-102
 Bone Holding Forceps
 185 mm, 7?",



Bone Rongeur
JU-54-103
 Bone Holding Forceps
 200 mm, 8",



Bone Rongeur
JU-54-104
 Bone Holding Forceps
 215 mm, 8?",



Bone Rongeur
JU-54-105
 Bone Rongeurs
 145 mm, 5?",



Friedman-Micro 6"
JU-54-106
 Bone Rongeurs
 150 mm, 6",

Bone Rongeurs


Friedman-Micro 5"
JU-54-107
 Bone Rongeurs
 145 mm, 5\ \ \ \,



Bone Rongeur
JU-54-108
 Bone Rongeurs
 140 mm, 5?",



LUER 6"
JU-54-109
 Bone Rongeurs
 150 mm, 6,



LUER 6"
JU-54-110
 Bone Rongeurs
 155 mm, 6,



2A-MEAD, SQUARE TIPS 5"
JU-54-111
 Bone Rongeur
 145 mm, 5,



BLUMENTHAL 6"
JU-54-112
 Bone Rongeur
 155 mm, 6,

Bone Rongeurs


BLUMENTHAL 45° 6"
JU-54-113
 Bone Rongeur
 150 mm, 6,



Bone Rongeur 7"
JU-54-114
 Bone Rongeur
 175 mm, 7",



Bone Rongeur
JU-54-115
 Bone Rongeur
 170 mm, 6 2/3",



BEYER 7"
JU-54-116
 Bone Rongeur
 180 mm, 7",



MEAD 6 1/2"
JU-54-117
 Bone Rongeur
 165 mm, 6 1/2",



BANE-HEARTMANN 7"
JU-54-118
 Bone Rongeur
 175 mm, 7",

Bone Rongeurs

BEYER DOUBLE-ACTION 7"
JU-54-119

 Bone Rongeur
 180 mm, 7",


ZUFAL-JANSEN 7"
JU-54-120

 Bone Rongeur
 180 mm, 7",

Bone Rongeur
JU-54-121

 Bone Rongeur
 130 mm, 5½",

Bone Rongeur
JU-54-122

 Bone Rongeur
 130 mm, 5½",

Cleveland 5"
JU-54-123

 Bone Cutting Forceps
 130 mm, 5",

Bone Rongeur
JU-54-124

 Bone Cutting Forceps
 150 mm, 6",

Bone Files




Fig. 1
Polokoff

Polokoff
JU-55-101
Bone Files
Fig- 1,




Fig. 33

JU-55-102
Bone Files
Fig- 33,



Fig. 64
Miler

Miler
JU-55-103
Bone Files
Fig- 64,




Fig. 45
Miler

Miler
JU-55-104
Bone Files
Fig- 45,



Fig. 10

JU-55-105
Bone Files
Fig- 10,



Fig. 1
Wahl

Wahl
JU-55-106
Bone Files
Fig- 1,

Bone Files



Fig. 2
Wahl

Wahl
JU-55-107
Bone Files
Fig- 2,



Fig. 12CA

JU-55-108
Bone Files
Fig- 12CA,



Fig. 12A

JU-55-109
Bone Files
Fig- 12A,




Fig. 12
Howard

Howard
JU-55-110
Bone Files
Fig- 12,

Bone Curettes



Fig 1
Hemingway

Hemingway
JU-56-101
Bone Curettes
Fig- 1,



Fig. 2
Hemingway

Hemingway
JU-56-102
Bone Curettes
Fig- 2,



Fig 3
Hemingway

Hemingway
JU-56-103
Bone Curettes
Fig- 3,



Fig 4
Hemingway

Hemingway
JU-56-104
Bone Curettes
Fig- 4,



Fig. 85
Lucas

Lucas
JU-56-105
Bone Curettes
Fig- 85,



Fig. 86
Lucas

Lucas
JU-56-106
Bone Curettes
Fig- 86,

Bone Curettes



Fig. 87
Lucas

Lucas
JU-56-107
Bone Curettes
Fig- 87,



Fig. 88
Lucas

Lucas
JU-56-108
Bone Curettes
Fig- 88,



Fig. 1
Hemingway

Hemingway
JU-56-109
Bone Curettes
Fig- 1,



Fig. 2
Hemingway

Hemingway
JU-56-110
Bone Curettes
Fig- 2,



Fig. 3
Hemingway

Hemingway
JU-56-111
Bone Curettes
Fig- 3,



Fig. 4
Hemingway

Hemingway
JU-56-112
Bone Curettes
Fig- 4,

Bone Curettes



Fig. 85
Lucas

Lucas
JU-56-113
Bone Curettes
Fig- 85,



Fig. 86
Lucas

Lucas
JU-56-114
Bone Curettes
Fig- 86,



Fig. 87
Lucas

Lucas
JU-56-115
Bone Curettes
Fig- 87,



Fig. 88
Lucas

Lucas
JU-56-116
Bone Curettes
Fig- 88,

Cotton & Dressing Tweezers

Cotton & Dressing Tweezers
JU-57-101

 Tooth Tweezers
 130 mm, 5 $\frac{1}{8}$ "

College Tweezer
JU-57-102

 Cotton & Dressing Tweezers
 130 mm, 5 $\frac{1}{8}$ "

College Tweezer Smooth Plan Handle
JU-57-103

 Cotton & Dressing Tweezers
 150 mm, 6"

Meriam Tweezers
JU-57-104

 Cotton & Dressing Tweezers
 150 mm, 6"

College Tweezer
JU-57-105

 Cotton & Dressing Tweezers
 160 mm, 6"

Meriam Tweezers
JU-57-106

 Cotton & Dressing Tweezers
 160 mm, 6"

Cotton & Dressing Tweezers



Meriam Tweezers
JU-57-107
 Cottn & Dressing Tweezers
 160 mm,6,



College Tweezer
JU-57-108
 Cottn & Dressing Tweezers
 150 mm,6,



College Tweezer
JU-57-109
 Cottn & Dressing Tweezers
 150 mm,6,



Self Locking Tweezer Serr.
JU-57-110
 Cottn & Dressing Tweezers
 150 mm,6,



College Tweezer
JU-57-111
 Cottn & Dressing Tweezers
 150 mm,6,



College Tweezer
JU-57-112
 Cottn & Dressing Tweezers
 150 mm,6,

Cotton & Dressing Tweezers



Perry Tweezer
JU-57-113

Cotton & Dressing Tweezers
130 mm,



JU-57-114

Cotton & Dressing Tweezers
150 mm,6,



Flagg Tweezer Smooth
JU-57-115

Cotton & Dressing Tweezers
150 mm,6,



Flagg Tweezer Serrated
JU-57-116

Cotton & Dressing Tweezers
150 mm,6,



JU-57-117

Cotton & Dressing Tweezers
150 mm,6,



Cotton & Dressing Tweezers
JU-57-118

Tooth Tweezers
150 mm,6",

Cotton & Dressing Tweezers



Cotton & Dressing Tweezers

JU-57-119

Tooth Tweezers
150 mm, 6",



Self Locking Tweezer Smooth

JU-57-120

Cotton & Dressing Tweezers
150 mm, 6",



Suture Forcep

JU-57-121

Cotton & Dressing Tweezers
150 mm, 6",

Rubber Dam Punch & Clamps



Rubber Dam Punch Ivory
JU-58-101
Rubber-Dam Instruments
160 mm,6,



Anisworth
JU-58-102
Rubber-Dam Instruments
160 mm,6?,



Brewer
JU-58-103
Rubber-Dam Instruments
170 mm,6?,



Palmer
JU-58-104
Rubber-Dam Instruments
175 mm,7,

Orthodontic Pliers



BIRD BEAK PLIER
JU-59-101
 ORTHODONTIC PLIERS, ANGLE



CUTTING PLIER FOR ORTHODONTICS
JU-59-102
 ORTHODONTIC PLIERS NO.1



CUTTING PLIERS FOR ORTHODONTIC
JU-59-103
 ORTHODONTIC PLIERS NO.2



DE LA ROSA ORTHO PLIER
JU-59-104
 ORTHODONTIC GROOVED PLIERS, DE LA ROSA



DE LA ROSA ORTHO PLIERS
JU-59-105
 ORTHODONTIC NON GROOVED PLIERS, DE LA ROSA



DISTAL END CUTTER
JU-59-106
 DISTAL END CUTTER

Orthodontic Pliers



DISTAL END CUTTER
JU-59-107 ITC



HARD WIRE CUTTER TC
JU-59-108 3TC



LIGHT WIRE CUTTER
JU-59-109
MICRO-MINIATURE WIRE CUTTER



LIGHT WIRE CUTTER
JU-59-110 ITC



MC KELLOPS PLIERS FOR ORTHODON
JU-59-111
ORTHODONTIC PLIERS, MC KELLOPS



MICRO DISTAL END CUTTER
JU-59-112
MICRO DISTAL WIRE CUTTER

Orthodontic Pliers



NANCE PLIERS FOR ORTHODONTICS
JU-59-113
 ORTHODONTIC PLIERS NANCE



ORHTODONTIC PLIERS
JU-59-114
 ORTHODONTIC PLIERS TO FORM
 LINGUAL ARCHS



ORTHODONTIC PLIER
JU-59-115
 ORTHODONTIC PLIERS TO REMOVE
 BRACKETS



ORTHODONTIC PLIER
JU-59-116
 ORTHODONTIC ATTACHMENT REMOVER
 PLIERS



ORTHODONTIC PLIERS
JU-59-117
 ORTHODONTIC ATTACHMENTS
 REMOVER PLIERS



ORTHODONTIC PLIERS
JU-59-118
 ORTHODONTIC DIRECT BONDING PLIERS

Orthodontic Pliers


PLIERS FOR ORTHODONTIC CURVED
JU-59-119
 ORTHODONTIC PLIERS HOW NO.2
 CURVED



PLIERS FOR ORTHODONTICS
JU-59-120
 ORTHODONTIC PLIERS



PLIERS FOR ORTHODONTICS
JU-59-121
 ORTHODONTIC PLIERS NO.2



PLIERS FOR ORTHODONTICS
JU-59-122
 ORTHODONTIC PLIERS WITH SMOOTH
 BEAKS



PLIERS FOR ORTHODONTICS
JU-59-123
 ORTHODONTIC PLIERS NO.2



PLIERS FOR ORTHODONTICS
JU-59-124
 ORTHODONTIC CONTOURING PLIERS,
 REYNOLDS

Orthodontic Pliers



PLIERS FOR ORTHODONTICS
JU-59-125
ORTHODONTIC PLIERS, WALDSACHS



PLIERS FOR ORTHODONTICS
JU-59-126
ORTHODONTIC PLIERS JOHNSON



PLIERS FOR ORTHODONTICS
JU-59-127
ORTHODONTIC PLIERS, NANCE



PLIERS FOR ORTHODONTICS
JU-59-128
ORTHODONTIC PLIERS, YOUNG



PLIERS FOR ORTHODONTICS
JU-59-129
ORTHODONTIC PLIERS



PLIERS FOR ORTHODONTICS
JU-59-130
ORTHODONTIC PLIERS WITH ROUNDED
JAW

Orthodontic Pliers



PLIERS FOR ORTHODONTICS
JU-59-131
ORTHODONTIC PLIERS, ANGLE



PLIERS FOR ORTHODONTICS
JU-59-132
ORTHODONTIC PLIERS, ADAMS



PLIERS FOR ORTHODONTICS
JU-59-133
ORTHODONTIC PLIERS WITH ROUND NOSE



PLIERS FOR ORTHODONTICS
JU-59-134
ORTHODONTIC PLIERS, ADERER



PLIERS FOR ORTHODONTICS
JU-59-135
ORTHODONTIC PLIERS ADERER



PLIERS FOR ORTHODONTICS
JU-59-136
ORTHODONTIC PLIERS, TWEED

Orthodontic Pliers



PLIERS FOR ORTHODONTICS
JU-59-137
 ORTHODONTIC PLIERS



PLIERS FOR ORTHODONTICS
JU-59-138
 ORTHODONTIC PLIERS, ANGLE NO.1



PLIERS FOR ORTHODONTICS STR.
JU-59-139
 ORTHODONTIC PLIERS HOW NO.1
 STRAIGHT



POSTERIOR BAND REMOVER
JU-59-140
 ORTHODONTIC BAND REMOVER PLIERS



POSTERIOR BAND REMOVER
JU-59-141
 ORTHODONTIC BAND REMOVER PLIERS



UNIVERSAL PLIERS FOR ORTHODON.
JU-59-142
 UNIVERSAL ORTHODONTIC PLIERS

Orthodontic Pliers



UNIVERSAL ORTHO PLIERS
JU-59-143
 UNIVERSAL BENDING ORTHODONTIC
 PLIERS



ROUND AND CONCAVE PLIERS
JU-59-144
 ORTHODONTIC PLIERS, ROUND AND
 CONCAVE



TWEED ORTHO PLIER
JU-59-145
 ORTHODONTIC PLIERS, TWEED



TWEED ORTHO PLIERS
JU-59-146
 ORTHODONTIC PLIERS, TWEED



TWEED ORTHO PLIERS
JU-59-147
 ORTHODONTIC LOOP FORMING
 PLIERS, TWEED



WEINGART ORTHO PLIERS
JU-59-148
 ORTHODONTIC UNIVERSAL PLIERS,
 WEINGART

Orthodontic Pliers



WIRE CUTTER

JU-59-151

WIRE CUTTER FOR ANY SIZE LIGATURE
WIRES



WIRE CUTTER

JU-59-152 ITC



HARD WIRE CUTTER TC

JU-59-153-3TC

Needle Holders


Crile-Murray
JU-60-101
 Needle Holders
 130 mm, 5,



Mayo-Hegar
JU-60-102
 Needle Holders
 145 mm, 5,



Crile-Murray
JU-60-103
 Needle Holders
 150 mm, 6,



Mayo-Hegar
JU-60-104
 Needle Holders
 150 mm, 6,




Mayo-Hegar
JU-60-105
 Needle Holders
 180 mm, 7,



Mayo-Hegar
JU-60-106
 Needle Holders
 200 mm, 8,

Needle Holders


Olsen-Hegar
JU-60-107
 Needle Holders
 140 mm,5,



Gillies
JU-60-108
 Needle Holders
 150 mm,6,



Boynton
JU-60-109
 Needle Holders
 125 mm,5,



Mathieu
JU-60-110
 Needle Holders
 170 mm,6,



Mathieu
JU-60-111
 Needle Holders
 200 mm,8,



Crile
JU-60-112
 Needle Holders
 150 mm,6,

Needle Holders


Toennis
JU-60-113
 Needle Holders
 180 mm, 7,



Mayo-Hegar
JU-60-114
 Needle Holders
 165 mm, 6,



Mayo-Hegar
JU-60-115
 Needle Holders
 210 mm, 8,



Mayo-Hegar
JU-60-116
 Needle Holders
 165 mm, 6,



Mayo-Hegar
JU-60-117
 Needle Holders
 210 mm, 8,



Debaquey
JU-60-118
 Needle Holders
 150 mm, 6,

Needle Holders


Debakey
JU-60-119
 Needle Holders
 180 mm, 7,



Ryder
JU-60-120
 Needle Holders
 135 mm, 5,



Ryder
JU-75-121
 Needle Holders
 155 mm, 6,



Ryder 7 CARB-N-SERT Jaws
JU-60-122
 Needle Holders
 175 mm, 7,



Bozemann 5
JU-60-123
 Needle Holders
 135 mm, 5,



Olsen-Hegar CARB-N-SERT Jaws
JU-60-124
 Needle Holders
 125 mm, 5,

Needle Holders


Olsen-Hegar CARB-N-SERT Jaws
JU-60-125
 Needle Holders
 165 mm, 6,



Olsen-Hegar CARB-N-SERT Jaws
JU-60-126
 Needle Holders
 160 mm, 6,



Olsen-Hegar CARB-N-SERT Jaws
JU-60-127
 Needle Holders
 145 mm, 5,



Needle Holders
JU-60-128
 Needle Holders
 185 mm, 7?"



Needle Holders
JU-60-129
 Needle Holders
 180 mm, 7"



Ryder 8 CARB-N-SERT Jaws
JU-60-130
 Needle Holders
 200 mm, 8,

Needle Holders


Needle Holders
JU-60-131
 Needle Holders
 200 mm, 8",



Mathieu 5
JU-60-132
 Needle Holders
 140 mm, 5,



Mathieu 6
JU-60-133
 Needle Holders
 170 mm, 6,



Mathieu 8
JU-60-134
 Needle Holders
 200 mm, 8,



Mathieu 5
JU-60-135
 Needle Holders
 140 mm, 5,



Mathieu 6
JU-60-136
 Needle Holders
 170 mm, 6,

Needle Holders



Mathieu 8
JU-60-137
Needle Holders
200 mm, 8,



Mathieu 7
JU-60-138
Needle Holders
195 mm, 7,



Mathieu CARB-N-SERT Jaws
JU-60-139
Needle Holders
205 mm, 8,

Gum Scissors



LISTER Bandage Scissors

JU-61-101
GUM Scissors
Extra fine
4-1/2" (11.4 cm),



LISTER Bandage Scissors

JU-61-102
Gum Scissors
Extra fine
5-1/2" (14 cm),



LISTER Bandage Scissors

JU-61-103
Gum Scissors
Extra fine
7-1/4" (18.4 cm),



KNOWLES Bandage Scissors

JU-61-104
Gum Scissors
Straight
5-1/2" (14 cm),



HEATH Suture Scissors

JU-61-105
Gum Scissors
Serrated blade
6-1/8" (15.6 cm),



SPENCER Stitch Scissors

JU-61-106
Gum Scissors
Delicate
3-1/2" (8.9 cm),

Gum Scissors



SHORTBENT Stitch Scissors

JU-61-107

Gum Scissors
Curved delicate
3-1/2" (8.9 cm),



LITTAUER JR. Stitch Scissors

JU-61-108

Gum Scissors
4-1/2" (11.4 cm),



NORTHBENT Stitch Scissors

JU-61-109

Gum Scissors
Curved light pattern
4-3/4" (12.1 cm),



LITTAUER Stitch Scissors

JU-61-110

Gum Scissors
Standard pattern
5-1/2" (14 cm),



LONG Oral Surgery Stitch Scissors

JU-61-111

Gum Scissors
6" (15.2 cm),



MILTEX Stitch Scissors

JU-61-112

Gum Scissors
Angled delicate
4-1/2" (11.4 cm),

Gum Scissors

REEH Stitch Scissors
JU-61-113

Gum Scissors

 Sharp tips small hook on 1 blade ribbon type
 3-3/4" (9.5 cm),

MILTEX Stitch Scissors
JU-61-114

Gum Scissors

 Curved sharp pointed tips ribbon type
 3-7/8" (9.8 cm),

MILTEX Stitch Scissors
JU-61-115

Gum Scissors

 Curved sharp pointed tips ribbon type
 4-1/2" (11.4 cm),

VANTAGE Crown Scissors
JU-61-116

Gum Scissors

 Collar curved
 4" (10.2 cm),

VANTAGE Collar
JU-61-117

Gum Scissors

 Crown Scissors straight
 4" (10.2 cm),

MILTEX Wire Cutting Scissors
JU-61-118

Gum Scissors

 Angled to side one serrated blade
 4-3/4" (12.1 cm),

Gum Scissors

MILTEX Wire Cutting Scissors
JU-61-119

Gum Scissors

 Angled to side, one serrated blade, Carb-N-Sert inserts
 4-3/4" (12.1 cm),

Standard Pattern Operating Scissors
JU-61-120

Gum Scissors

 Straight sharp-sharp points
 4-1/2" (11.4 cm),

Standard Pattern Operating Scissors
JU-61-121

Gum Scissors

 Straight sharp-sharp points
 5-1/2" (14 cm),

Standard Pattern Operating Scissors
JU-61-122

Gum Scissors

 Straight sharp-blunt points
 4-1/2" (11.4 cm),

Standard Pattern Operating Scissors
JU-61-123

Gum Scissors

 Straight sharp-blunt points
 5-1/2" (14 cm),

Standard Pattern Operating Scissors
JU-61-124

Gum Scissors

 Straight sharp-blunt points
 5-1/2" (14 cm),

Gum Scissors

Standard Pattern Operating Scissors
JU-61-125

Gum Scissors
Curved sharp/sharp points
4-1/2" (11.4 cm),


Standard Pattern Operating Scissors
JU-61-126

Gum Scissors
Curved sharp/sharp points
5-1/2" (14 cm),


Standard Pattern Operating Scissors
JU-61-127

Gum Scissors
Curved
4-1/2,


Standard Pattern Operating Scissors
JU-61-128

Gum Scissors
Curved, sharp-blunt points
5-1/2" (14 cm),


MAYO Dissecting Scissors
JU-61-129

Gum Scissors
Straight, standard beveled blades
5-1/2" (14 cm),


MAYO Dissecting Scissors
JU-61-130

Gum Scissors
Straight, standard beveled blades
5-1/2" (14 cm),

Gum Scissors

MAYO Dissecting Scissors
JU-61-131

Gum Scissors
Curved, standard beveled blades
5-1/2" (14 cm),


MAYO Dissecting Scissors
JU-61-132

Gum Scissors
Curved, standard beveled blades
5-1/2" (14 cm),


Delicate Pattern METZENBAUM Scissors
JU-61-133

Gum Scissors
Straight, blunt points
5-1/2" (14 cm),


Delicate Pattern METZENBAUM Scissors
JU-61-134

Gum Scissors
Straight, sharp points
5-1/2" (14 cm),


Standard Pattern METZENBAUM Scissors
JU-61-135

Gum Scissors
Curved, blunt points
5-1/2" (14 cm),


Standard Pattern METZENBAUM Scissors
JU-61-136

Gum Scissors
Curved, blunt points
5-1/2" (14 cm),

Gum Scissors

Delicate Pattern METZENBAUM Scissors
JU-61-137

Gum Scissors
Curved, blunt points
5-1/2" (14 cm),


Delicate Pattern METZENBAUM Scissors
JU-61-138

Gum Scissors
Curved, sharp points
5-1/2" (14 cm),


KELLY Scissors
JU-61-139

Gum Scissors
Straight, sharp points
6-1/4" (15.9 cm),


KELLY Scissors
JU-61-140

Gum Scissors
Curved, sharp points
6-1/4" (15.9 cm),


KELLY Scissors
JU-61-141

Gum Scissors
Sharp points, angled on side
6-1/4" (15.9 cm),


KAYE Blepharoplasty and Dissecting Scissors
JU-61-142

Gum Scissors
One serrated blade, ribbon type
4-1/4" (10.8 cm),

Gum Scissors



DEAN Dissecting Scissors

JU-61-143

Gum Scissors

Blades angled on flat, one serrated blade
6-3/4" (17.1 cm),



KAYE Dissecting Scissors

JU-61-144

Gum Scissors

Curved, one serrated blade
4-17/32" (11.5 cm),



WAGNER Plastic Surgery Scissors

JU-61-145

Gum Scissors

Straight, sharp-sharp points, serrated blade
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

JU-61-146

Gum Scissors

Straight, sharp-sharp points, serrated blade
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

JU-61-147

Gum Scissors

Straight, sharp-blunt points
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

JU-61-148

Gum Scissors

Straight, sharp-blunt points, serrated blade
4-3/4" (12.1 cm),

Gum Scissors

WAGNER Plastic Surgery Scissors
JU-61-149

Gum Scissors
Straight, blunt-blunt points
4-3/4" (12.1 cm),


WAGNER Plastic Surgery Scissors
JU-61-150

Gum Scissors
Straight, blunt-blunt points, serrated blade
4-3/4" (12.1 cm),


WAGNER Plastic Surgery Scissors
JU-61-151

Gum Scissors
Curved, sharp-sharp points
4-3/4" (12.1 cm),


WAGNER Plastic Surgery Scissors
JU-61-152

Gum Scissors
Curved, sharp-sharp points, serrated blade
4-3/4" (12.1 cm),


WAGNER Plastic Surgery Scissors
JU-61-153

Gum Scissors
Curved, sharp-blunt points
4-3/4" (12.1 cm),


WAGNER Plastic Surgery Scissors
JU-61-154

Gum Scissors
Curved, sharp-blunt points, serrated blade
4-3/4" (12.1 cm),

Gum Scissors

WAGNER Plastic Surgery Scissors
JU-61-155

Gum Scissors
Curved, blunt-blunt points
4-3/4" (12.1 cm),


WAGNER Plastic Surgery Scissors
JU-61-156

Gum Scissors
Curved, blunt-blunt points, serrated blade
4-3/4" (12.1 cm),


WAGNER Plastic Surgery Scissors
JU-61-157

Gum Scissors
Straight, sharp-sharp points
4-3/4" (12.1 cm),


WAGNER Plastic Surgery Scissors
JU-61-158

Gum Scissors
Straight, sharp-sharp point, serrated blade
4-3/4" (12.1 cm),


Delicate Pattern METZENBAUM Scissors
JU-61-159

Gum Scissors
Curved
5" (12.7 cm),


Iris Scissors
JU-61-160

Gum Scissors
Straight
4" (10.2 cm),

Gum Scissors


Iris Scissors
JU-61-161
 Gum Scissors
 Curved
 4" (10.2 cm),



Iris Scissors
JU-61-162
 Gum Scissors
 Straight
 4-1/2" (11.4 cm),



Iris Scissors
JU-61-163
 Gum Scissors
 Straight
 4-1/2" (11.4 cm),



Iris Scissors
JU-61-164
 Gum Scissors
 Curved
 4-1/2" (11.4 cm),



Iris Scissors
JU-61-165
 Gum Scissors
 Curved
 4-1/2" (11.4 cm),



Iris Scissors
JU-61-166
 Gum Scissors
 Angled on side
 4-1/2" (11.4 cm),

Gum Scissors

BROPHY Scissors
JU-61-167

Gum Scissors
Straight, sharp points
5-1/2" (14 cm),


BROPHY Scissors
JU-61-168

Gum Scissors
Curved, sharp points
5-1/2" (14 cm),


GOLDMAN-FOX Scissors
JU-61-169

Gum Scissors
Curved fine tips, one serrated blade
4" (10.2 cm),


GOLDMAN-FOX Scissors
JU-61-170

Gum Scissors
Curved fine tips, one serrated blade
5-1/2" (14 cm),


KELLY Scissors
JU-61-171

Gum Scissors
Sharp points, straight, one serrated blade
6-1/4" (15.9 cm),


KELLY Scissors
JU-61-172

Gum Scissors
Sharp points, curved, one serrated blade
6-1/4" (15.9 cm),

Gum Scissors

KELLY Scissors
JU-61-173

Gum Scissors

 Sharp points, angular, one serrated blade
 6-1/4" (15.9 cm),

Iris Operating Scissors
JU-61-174

Gum Scissors

 Curved sideways
 4-1/2" (11.4 cm),

GOLDMAN-FOX Scissors
JU-61-175

Gum Scissors

 Curved, one serrated blade, angled shank
 4-1/2" (11.4 cm),

LA GRANGE Scissors
JU-61-176

Gum Scissors

 Curved, one serrated blade, angled shank
 4-1/4" (10.8 cm),

BROPHY Operating Scissors
JU-61-177

Gum Scissors

 Curved, sharp/sharp, with one serrated blade
 5-1/2" (14 cm),

LOCKLIN Operating Scissors
JU-61-178

Gum Scissors

 Angled on side 25 degrees, one serrated
 blade
 6-1/4" (15.9 cm),

Gum Scissors

LOCKLIN Operating Scissors
JU-61-179

Gum Scissors

 Curved shank, angled on side 25 degrees,
 one serrated blade
 6-1/4" (15.9 cm),

MILTEX SuperCut GOLDMAN-FOX Gum Scissors
JU-61-180

Gum Scissors

 Curved
 4-1/2" (11.4 cm),

HEATH Suture Scissors
JU-61-181

Gum Scissors

 Serrated blade
 6-1/8" (15.6 cm),

Iris Scissors
JU-61-182

Gum Scissors

 Straight, with 20 mm blades, delicate
 3-1/2" (8.9 cm),

Iris Scissors
JU-61-183

Gum Scissors

 Curved, with 20 mm blades, delicate
 3-1/2" (8.9 cm),

NOYES (ROSS) Iris Scissors
JU-61-184

Gum Scissors

 Curved, sharp/sharp points
 4-1/2" (11.4 cm),

Gum Scissors

LOCKLIN Operating Scissors
JU-61-179

Gum Scissors

 Curved shank, angled on side 25 degrees,
 one serrated blade
 6-1/4" (15.9 cm),

MILTEX SuperCut GOLDMAN-FOX Gum Scissors
JU-61-180

Gum Scissors

 Curved
 4-1/2" (11.4 cm),

HEATH Suture Scissors
JU-61-181

Gum Scissors

 Serrated blade
 6-1/8" (15.6 cm),

Iris Scissors
JU-61-182

Gum Scissors

 Straight, with 20 mm blades, delicate
 3-1/2" (8.9 cm),

Iris Scissors
JU-61-183

Gum Scissors

 Curved, with 20 mm blades, delicate
 3-1/2" (8.9 cm),

NOYES (ROSS) Iris Scissors
JU-61-184

Gum Scissors

 Curved, sharp/sharp points
 4-1/2" (11.4 cm),

Gum Scissors

CASTROVIEJO Corneal Scissors
JU-61-185

Gum Scissors
Curved, sharp points
3-3/4" (9.5 cm),


QUIMBY Gum Scissors
JU-61-186

Gum Scissors
Strongly curved and pointed, smooth blades
15 mm
4" (10.2 cm),


JOSEPH Scissors
JU-61-187

Gum Scissors
Straight
5-3/4" (14.6 cm),


JOSEPH Scissors
JU-61-188

Gum Scissors
Curved
5-3/4" (14.6 cm),

Amalgam Instruments



Fig. SHO-A

JU-62-117
Amalgam Instruments
Fig- SHO-A,



Fig. 3-SHO

JU-62-118
Amalgam Instruments
Fig- 3-SHO,



Fig. 104

JU-62-119
Amalgam Instruments
Fig- 104,



Fig. 8

JU-62-120
Amalgam Instruments
Fig- 8,



Fig. 1

JU-62-121
Amalgam Instruments
Fig- 1,



Fig. 2

JU-62-122
Amalgam Instruments
Fig- 2,

Amalgam Instruments



Fig 3

JU-62-123
Amalgam Instruments
Fig- 3,



Fig 4

JU-62-124
Amalgam Instruments
Fig- 4,



Fig. SHO-A

JU-62-125
Amalgam Instruments
Fig- SHO-A,



Fig. 3-SHO

JU-62-126
Amalgam Instruments
Fig- 3-SHO,



Fig. 104

JU-62-127
Amalgam Instruments
Fig- 104,



Fig. 8

JU-62-128
Amalgam Instruments
Fig. 8,

Amalgam Instruments



JU-62-129

Amalgam Instruments
Fig- 1,



JU-62-130

Amalgam Instruments
Fig- 2,



JU-62-131

Amalgam Instruments
Fig- 3,



JU-62-132

Amalgam Instruments
Fig- 4,



JU-62-133

Amalgam Instruments
Fig- 1,



JU-62-134

Amalgam Instruments
Fig- 2,

Amalgam Instruments



Fig. 3

JU-62-135

Amalgam Instruments
Fig- 3,



Fig. 4

JU-62-136

Amalgam Instruments
Fig- 4,



Ivory Fig. 8
5 mm
6 mm
7 mm

Toffemire-Universal

JU-62-137

Amalgam Instruments
Lvory Fig- 8,5 mm,6 mm,7 mm,



Nystrom I + II

Ivory Fig. 1

Ivory Fig. 12

JU-62-142

Amalgam Instruments
Lvory Fig- 1,Lvory Fig- 2,



Amalgam Instruments

JU-62-143

Amalgam Pistolen
210 mm,



Amalgam Instruments

JU-62-144

Amalgam Pistolen
185 mm,

Amalgam Instruments



Amalgam Instruments
JU-62-145
 Amalgam Pistolen
 135 mm,



Amalgam Instruments
JU-62-146
 Amalgam Pistolen
 145 mm,



Amalgam Carrier
JU-62-147
 ?
 Fig. 1, Fig. 2, Fig. 3,



Amalgam Carrier
JU-62-148
 ?
 Fig. 1, Fig. 2, Fig. 3,



Amalgam Carrier
JU-62-149
 ?
 Fig. 1, Fig. 2, Fig. 3,

Mouth Gags & Retractor



Heister
JU-63-141
Mouth Gags,
130 mm, 5 $\frac{1}{8}$,



MOLT PEDIATRIC
JU-63-142
Mouth Gags,
100 mm, 4,



MOLT CHILD
JU-63-143
Mouth Gags,
120 mm, 4,



MOLT ADULT
JU-63-144
Mouth Gags,
140 mm, 5,



FERGUSON-ACKLAND
JU-63-145
Mouth Gags,
145 mm, 5,



DENHART
JU-63-146
Mouth Gags,
125 mm, 5,

Mouth Gags & Retractor



ROSER-KONIG
JU-63-147
 Mouth Gags,
 160 mm, 6 1/4,




ROSER-KONIG
JU-63-148
 Mouth Gags,
 185 mm, 7,



LANGENBECK
JU-63-149
 Finger Protectors
 185 mm, 7,



JU-63-150
 Tongue Depressor and Cheek Retractors



Retractor
JU-63-151
 Mouth Gags,



HARTMANN
JU-63-152
 Retractors
 145 mm, 5,

Mouth Gags & Retractor



Retractor
JU-63-153
Mouth Gags,
140 mm,5,



SCHMID
JU-63-154
Mouth Gags,
140 mm,5,



Retractor
JU-63-155
Mouth Props
140 mm,5?,



Retractor
JU-63-156
Retractors
160 mm,6?,



VESTIBULUM RETRACTOR
JU-63-157
Retractors
140 mm,5,



STERNBERG
JU-63-158
Retractors
120 mm,4,

Mouth Gags & Retractor


ROUX
JU-63-159
 Retractors
 19x22 mm, 23x30 mm,



ROUX
JU-63-160
 Retractors
 25x26 mm, 28x38 mm,



ROUX
JU-63-161
 Retractors
 28x29 mm, 33x45 mm,



MIDDELDORPF
JU-63-162
 Retractors
 15x15, 200 mm, 8,



MIDDELDORPF
JU-63-163
 Retractors
 22x22 mm, 210 mm, 8,



MIDDELDORPF
JU-63-164
 Cheek Retractors
 26x30 mm, 230 mm, 9,

Mouth Gags & Retractor


Gauze Case
JU-63-165
 Gauze Case
 26x30 mm, 230 mm, 9,



Retractor
JU-63-166
 Retractors
 165 mm, 6?",




Retractor
JU-63-167
 Retractors
 165 mm, 6?",



Retractor
JU-63-168
 Retractors
 165 mm, 6?",



Retractor
JU-63-169
 Retractors
 165 mm, 6?",



Retractor
JU-63-170
 Retractors
 165 mm, 6?",

Mouth Gags & Retractor


**Retractor
JU-63-171**
Retractors
165 mm, 6 7/8",



**Retractor
JU-63-172**
Retractors
165 mm, 6 7/8",



**Retractor
JU-63-173**
Retractors
165 mm, 6 7/8",



**Retractor
JU-63-174**
Retractors
160 mm, 6 3/8",



**Retractor
JU-63-175**
Retractors
160 mm, 6 3/8",



**LANGENBECK-GREEN
JU-63-176**
Retractors
160 mm, 16x6 mm,

Periodontia Instruments



Fig. 1/2
105/100
Kirkland

**Kirkland
JU-64-101**

Periodontia Instruments
Fig -1/2,



Fig. 1/2
Orban

**Orban
JU-64-102**

Periodontia Instruments
Fig- 1/2,



Fig. 1/2
105/100
Kirkland

**Kirkland
JU-64-103**

Periodontia Instruments
Fig- 1/2,



Fig. 1/2
Orban

**Orban
JU-64-104**

Periodontia Instruments
Fig- 1/2,



Fig. 7
Goldman-Fox

**Goldman-Fox
JU-64-105**

Periodontia Instruments
Fig- 7,



Fig. 8
Goldman-Fox

**Goldman-Fox
JU-64-106**

Periodontia Instruments
Fig- 8,

Periodontia Instruments



Fig. 9
Goldman-Fox

Goldman-Fox
JU-64-107

Periodontia Instruments
Fig- 9,



Fig. 11
Goldman-Fox

Goldman-Fox
JU-64-108

Periodontia Instruments
Fig- 11,



Fig. 3/4
Buck

Buck
JU-64-109

Periodontia Instruments
Fig- 3/4,



Fig. 5/6
Buck

Buck
JU-64-110

Periodontia Instruments
Fig- 5/6,



Fig. 7
Goldman-Fox

Goldman-Fox
JU-64-111

Periodontia Instruments
Fig- 7,



Fig. 8
Goldman-Fox

Goldman-Fox
JU-64-112

Periodontia Instruments
Fig- 8,

Periodontia Instruments



Fig. 9
Goldman-Fox

Goldman-Fox
JU-64-113
Periodontia Instruments
Fig- 9,



Fig. 11
Goldman-Fox

Goldman-Fox
JU-64-114
Periodontia Instruments
Fig- 11,



Fig. 3/4
Buck

Buck
JU-64-115
Periodontia Instruments
Fig- 3/4,



Fig. 5/6
Buck

Buck
JU-64-116
Periodontia Instruments
Fig- 5/6,



Fig. 1/2
Sanders

Sanders
JU-64-117
Periodontia Instruments
Fig- 1/2,



Fig. 4/5
Sanders

Sanders
JU-64-118
Periodontia Instruments
Fig- 4/5,

Periodontia Instruments



Fig. 19/20

OCHSENBEIN CHISEL NO 4

JU-64-119

Periodontia Instruments
Fig- 19/20,



Fig. 5/6
Back

Crane-Kaplan

JU-64-120

Periodontia Instruments
Fig- CK3,



Fig. 1/2
Back

Sanders

JU-64-121

Periodontia Instruments
Fig- 1/2,



Fig. 4/5
SANDERS

Sanders

JU-64-122

Periodontia Instruments
Fig- 4/5,



Fig. 19/20

OCHSENBEIN CHISEL NO 4

JU-64-123

Periodontal Instruments
Fig- 19/20,



Fig. CK3
Crane-Kaplan

Crane-Kaplan

JU-64-124

Periodontia Instruments
Fig- CK3,

Periodontia Instruments



Fig. HU 1/2
Curette

Curette

JU-64-125

Periodontia Instruments
Fig- HU 1/2,



Fig. 2R/2L
Columbia

Columbia

JU-64-126

Periodontia Instruments
Fig- 2R/2L,



Fig. 4R/4L
Columbia

Columbia

JU-64-127

Periodontia Instruments
Fig- 4R/4L,



Fig. 13/14
Columbia

Columbia

JU-64-128

Periodontia Instruments
Fig- 13/14,



Fig. 411/412

JU-64-129

Periodontia Instruments
Fig- 411/412,



Fig. 413/414

JU-64-130

Periodontia Instruments
Fig- 413/414,

Periodontia Instruments



Zerling-Chiles

Zerling-Chiles
JU-64-131

Periodontia Instruments
Fig- 413/414,



Fig. HU 1/2
Curette

Curette
JU-64-132

Periodontia Instruments
Fig- HU 1/2,



Fig. 2R/2L
Columbia

Columbia
JU-64-133

Periodontia Instruments
Fig- 2R/2L,



Fig. 4R/4L
Columbia

Columbia
JU-64-134

Periodontia Instruments
Fig- 4R/4L,



Fig. 13/14
Columbia

Columbia
JU-64-135

Periodontia Instruments
Fig- 13/14,



Fig. 1/2
Langer

Langer
JU-64-136

Periodontia Instruments
Fig- 1/2,

Periodontia Instruments



Fig. 3/4
Langer

Langer

JU-64-137

Periodontia Instruments
Fig- 3/4,



Fig. 1/2
Barnhart

Barnhart

JU-64-138

Periodontia Instruments
Fig- 1/2,



Fig. 5/6
Barnhart

Barnhart

JU-64-139

Periodontia Instruments
Fig- 5/6,



Fig. ME-1
(15/11)

JU-64-140

Periodontia Instruments
Fig- ME-1,



Fig. 3/4
Langer (23/24)

JU-64-141

Periodontia Instruments
Fig- 3/4,



Fig. 1/2
Langer

Langer

JU-64-142

Periodontia Instruments
Fig- 1/2,

Periodontia Instruments



Fig. 3/4
Langer

Langer

JU-64-143

Periodontia Instruments
Fig- 3/4,



Fig. 1/2
Barnhart

Barnhart

JU-64-144

Periodontia Instruments
Fig- 1/2,



Fig. 5/8
Barnhart

Barnhart

JU-64-145

Periodontia Instruments
Fig- 5/6,



Fig ME-1
[20/11]

JU-64-146

Periodontia Instruments
Fig- ME-1,



Fig. 3/4
Rise [23/24]

JU-64-147

Periodontia Instruments
Fig- 3/4,



Fig. 1/2
Gracey

Gracey

JU-64-148

Periodontia Instruments
Fig- 1/2,

Periodontia Instruments



Fig. 3/4
Gracey

Gracey

JU-64-149

Periodontia Instruments
Fig- 3/4,



Fig. 5/6
Gracey

Gracey

JU-64-150

Periodontia Instruments
Fig- 5/6,



Fig. 7/8
Gracey

Gracey

JU-64-151

Periodontia Instruments
Fig- 7/8,



Fig. 9/10
Gracey

Gracey

JU-64-152

Periodontia Instruments
Fig- 9/10,



Fig. 11/12
Gracey

Gracey

JU-64-153

Periodontia Instruments
Fig- 11/12,



Fig. 13/14
Gracey

Gracey

JU-64-154

Periodontia Instruments
Fig- 13/14,

Periodontia Instruments



Fig. 1/2
Gracey

Gracey

JU-64-155

Periodontia Instruments
Fig- 1/2,



Fig. 3/4
Gracey

Gracey

JU-64-156

Periodontia Instruments
Fig- 3/4,



Fig. 5/6
Gracey

Gracey

JU-64-157

Periodontia Instruments
Fig- 5/6,



Fig. 7/8
Gracey

Gracey

JU-64-158

Periodontia Instruments
Fig- 7/8,



Fig. 9/10
Gracey

Gracey

JU-64-159

Periodontia Instruments
Fig- 9/10,



Fig. 11/12
Gracey

Gracey

JU-64-160

Periodontia Instruments
Fig- 11/12,

Periodontia Instruments



Fig. 13/14
Gracey

Gracey

JU-64-161

Periodontia Instruments
Fig- 13/14,



Fig. 11/14
Gracey

Gracey

JU-64-162

Periodontia Instruments
Fig- 11/14,



Fig. 12/13
Gracey

Gracey

JU-64-163

Periodontia Instruments
Fig- 12/13,



Fig. 1
De Marco

Ce Marco

JU-64-164

Periodontia Instruments
Fig- 1,



Fig. 2
De Marco

De Marco

JU-64-165

Periodontia Instruments
Fig- 2,



Fig. 13/14
Anterior

Anterior

JU-64-166

Periodontia Instruments
Fig- 13/14,

Periodontia Instruments



Fig. 17/18
Posterior

Posterior
JU-64-167

Periodontia Instruments
Fig- 17/18,



Fig. 11/14
Gracey

Gracey
JU-64-168

Periodontia Instruments
Fig- 11/14,



Fig. 12/13
Gracey

Gracey
JU-64-169

Periodontia Instruments
Fig- 12/13,



Fig. 1
De Marco

De Marco
JU-64-170

Periodontia Instruments
Fig- 1,



Fig. 2
De Marco

De Marco
JU-64-171

Periodontia Instruments
Fig- 2,



Fig. 13/14
Anterior

Anterior
JU-64-172

Periodontia Instruments
Fig- 13/14,

Periodontia Instruments



Fig. 17/18
Posierior

**Posierior
JU-64-173**

Periodontia Instruments
Fig- 17/18,



Fig. 1
Goldman-Fox

**Goldman-Fox
JU-64-174**

Periodontia Instruments
Fig- 1,



Fig. 2
Goldman-Fox

**Goldman-Fox
JU-64-175**

Periodontia Instruments
Fig- 2,



Fig. 3
Goldman-Fox

**Goldman-Fox
JU-64-176**

Periodontia Instruments
Fig- 3,



Fig. 4
Goldman-Fox

**Goldman-Fox
JU-64-177**

Periodontia Instruments
Fig- 4,



Fig. 5
Goldman-Fox

**Goldman-Fox
JU-64-178**

Periodontia Instruments
Fig- 5,

Periodontia Instruments



Fig. 6
Goldman-Fox

Goldman-Fox
JU-64-179
Periodontia Instruments
Fig- 6,



Fig. 21
Goldman-Fox

Goldman-Fox
JU-64-180
Periodontia Instruments
Fig- 21,



Fig. 1
Goldman-Fox

Goldman-Fox
JU-64-181
Periodontia Instruments
Fig- 1,



Fig. 2
Goldman-Fox

Goldman-Fox
JU-64-182
Periodontia Instruments
Fig- 2,



Fig. 3
Goldman-Fox

Goldman-Fox
JU-64-183
Periodontia Instruments
Fig- 3,



Fig. 4
Goldman-Fox

Goldman-Fox
JU-64-184
Periodontia Instruments
Fig- 4,

Periodontia Instruments



Fig. 5
Goldman-Fox

Goldman-Fox
JU-64-185
Periodontia Instruments
Fig- 5,



Fig. 6
Goldman-Fox

Goldman-Fox
JU-64-186
Periodontia Instruments
Fig- 6,



Fig. 21
Goldman-Fox

Goldman-Fox
JU-64-187
Periodontia Instruments
Fig- 21,



Fig. 11/12
Mc Call

Mc Call
JU-64-188
Periodontia Instruments
Fig- 11/12,



Fig. 11A/12A
Mc Call

Mc Call
JU-64-189
Periodontia Instruments
Fig- 11A/12A,



Fig. 13/14
Mc Call

Mc Call
JU-64-190
Periodontia Instruments
Fig- 13/14,

Periodontia Instruments


Fig. 13S/14S
Mc Call

Mc Call

JU-64-191

Periodontia Instruments
Fig- 13S/14S,



Fig. 17S/18S
Mc Call

Mc Call

JU-64-192

Periodontia Instruments
Fig- 17S/18S,



Fig. 19/20
Mc Call

Mc Call

JU-64-193

Periodontia Instruments
Fig- 19/20,



Fig. 11/12
Mc Call

Mc Call

JU-64-194

Periodontia Instruments
Fig- 11/12,



Fig. 11A/12A
Mc Call

Mc Call

JU-64-195

Periodontia Instruments
Fig- 11A/12A,



Fig. 13/14
Mc Call

Mc Call

JU-64-196

Periodontia Instruments
Fig- 13/14,

Periodontia Instruments



Fig. 13S/14S
Mc Call

Mc Call

JU-64-197

Periodontia Instruments
Fig- 13S/14S,



Fig. 17S/18S
Mc Call

Mc Call

JU-64-198

Periodontia Instruments
Fig- 17S/18S,



Fig. 19/20
Mc Call

Mc Call

JU-64-199

Periodontia Instruments
Fig- 19/20,



Fig. 2/3
Taylor

Taylor

JU-64-200

Periodontia Instruments
Fig- 2/3,



Fig. 204

JU-64-201

Periodontia Instruments
Fig- 204,



Fig. 204s

JU-64-202

Periodontia Instruments
Fig- 204s,

Periodontia Instruments



Fig. 204S

JU-64-203

Periodontia Instruments
Fig- 204S,



Fig. 6

Crane-Kaplan

Crane-Kaplan

JU-64-204

Periodontia Instruments
Fig- 6,



Fig. 7/8

Younger-Good

Younger-Good

JU-64-205

Periodontia Instruments
Fig- 7/8,



Fig. 2/3

Taylor

Taylor

JU-64-206

Periodontia Instruments
Fig- 2/3,



Fig. 204

JU-64-207

Periodontia Instruments
Fig- 204,



Fig. 204S

JU-64-208

Periodontia Instruments
Fig- 204S,

Periodontia Instruments



Fig. 204S

JU-64-209
Periodontia Instruments
Fig- 204S,



Fig. 6
Crane-Kaplan

Crane-Kaplan
JU-64-210
Periodontia Instruments
Fig- 6,



Fig. 7/8
Younger-Good

Younger-Good
JU-64-211
Periodontia Instruments
Fig- 7/8,



Fig. 3/7
Hirschfeld

Hirschfeld
JU-64-212
Periodontia Instruments
Fig- 3/7,



Fig. 30/33
Jaquette

Jaquette
JU-64-213
Periodontia Instruments
Fig- 30/33,



Fig. 31/32
Jaquette

Jaquette
JU-64-214
Periodontia Instruments
Fig- 31/32,

Periodontia Instruments



Fig. 31/32
Jaquette

Jaquette
JU-64-215

Periodontia Instruments
Fig- 31/32,



Fig. 3/7
Hirschfeld

Hirschfeld
JU-64-216

Periodontia Instruments
Fig- 3/7,



Fig. 30/33
Jaquette

Jaquette
JU-64-217

Periodontia Instruments
Fig- 30/33,



Fig. 31/22
Jaquette

Jaquette
JU-64-218

Periodontia Instruments
Fig- 31/32,



Fig. 34/35
Jaquette

Jaquette
JU-64-219

Periodontia Instruments
Fig- 34/35,



Fig. 37/38
Jaquette

Jaquette
JU-64-220

Periodontia Instruments
Fig- 37/38,

Periodontia Instruments



Fig. 1
Ochsenbein

Ochsenbein
JU-64-221
Periodontia Instruments
Fig- 1,



Fig. 2
Ochsenbein

Ochsenbein
JU-64-222
Periodontia Instruments
Fig- 2,



Fig. 1
Gardner

Gardner
JU-64-223
Periodontia Instruments
Fig- 1,



Fig. 2
Gardner

Gardner
JU-64-224
Periodontia Instruments
Fig- 2,



Fig. 5
Gardner

Gardner
JU-64-225
Periodontia Instruments
Fig- 5,



Fig. 5R
Gardner

Gardner
JU-64-226
Periodontia Instruments
Fig- 5R,

Periodontia Instruments



Fig. 6
Gardner

Gardner
JU-64-227

Periodontia Instruments
Fig- 6,



Fig. 1W
Wakef Old

Wakef Old
JU-64-228

Periodontia Instruments
Fig- 1W,



Fig. 2W
Wakef Old

Wakef Old
JU-64-229

Periodontia Instruments
Fig- 2W,



Fig. 3
Buckley

Buckley
JU-64-230

Periodontia Instruments
Fig- 3,



Fig. M9
Molt

Molt
JU-64-231

Periodontia Instruments
Fig- M9,



Fig. PR3
Prichard

Prichard
JU-64-232

Periodontia Instruments
Fig- PR3,

Periodontia Instruments



Fig. 1
Hopkins

Hopkins
JU-64-233

Periodontia Instruments
Fig- 1,



Fig. 2
Hopkins

Hopkins
JU-64-234

Periodontia Instruments
Fig- 2,



Fig. 15
Frerr

Frerr
JU-64-235

Periodontia Instruments
Fig- 15,



Fig. 16
Frerr

Frerr
JU-64-236

Periodontia Instruments
Fig- 16,



Fig. 3
Mead

Mead
JU-64-237

Periodontia Instruments
Fig- 3,



Fig. 23
Seldin retractor

Seldin retractor
JU-64-238

Periodontia Instruments
Fig- 23,

Periodontia Instruments



Williger
cm. 17

**Williger
JU-64-239**

Periodontia Instruments
cm 17,



Williger
cm. 16

**Williger
JU-64-240**

Periodontia Instruments
cm 16,



Williger
cm. 16

**Williger
JU-64-241**

Periodontia Instruments
cm 16,



Williger
cm. 13

**Williger
JU-64-242**

Periodontia Instruments
cm 13,



Fig. 1S/2S
Sugarman

**Sugarman
JU-64-243**

Periodontia Instruments
Fig- 1S/2S,



Fig. 3S/4S
Sugarman

**Sugarman
JU-64-244**

Periodontia Instruments
Fig- 3S/4S,

Periodontia Instruments



Fig. 1/2
TG Chisel

TG Chsel
JU-64-245

Periodontia Instruments
Fig- 1/2,



Fig. 1/2
TG Chisel

TG Chsel
JU-64-246

Periodontia Instruments
Fig- 1/2,



Fig. 1
TG Chisel

TG Chsel
JU-64-247

Periodontia Instruments
Fig- 1,



Fig. 2
TG Chisel

TG Chsel
JU-64-248

Periodontia Instruments
Fig- 2,



Fig. 1
(15K)
Kirkland

Kirkland
JU-64-249

Periodontia Instruments
Fig- 1,



Fig. 2
(16K)
Kirkland

Kirkland
JU-64-250

Periodontia Instruments
Fig- 2,

Periodontia Instruments



Fig. 1
Orban

Orban
JU-64-251
Periodontia Instruments
Fig- 1,



Fig. 2
Orban


Orban
JU-64-252
Periodontia Instruments
Fig- 2,

Filling instruments




Westcott
JU-65-101

Filling instruments
Fig- 1,




Westcott
JU-65-110

Filling instruments
Fig- 1,



JU-65-119

Filling instruments
Fig- 1, Fig- 2, Fig- 3, Fig- 4, Fig- 5,



JU-65-125

Filling instruments
Fig- 1, Fig- 2, Fig- 3, Fig- 1L, Fig- 2L, Fig- 3L,

Cement Spatulas



JU-66-101
No. 22 Cement Spatula



JU-66-102
No. 24A Cement Spatula, Rigid



JU-66-103
No. 24 Cement Spatula, Flexible



JU-66-104
No. 313 Cement Spatula,




JU-66-105
No. 5 Cement Spatula




JU-66-106
No. 7 Cement Spatula

Single/Double End Explorers



JU-67-101

Single Ended/ Double Ended Explorers
Fig- UT 1/2,



JU-67-112

Single Ended/ Double Ended Explorers
Fig- UT 1/2,

Mirror Handles



SIMPLE STEM PLANE 4
JU-68-113
MOUTH MIRROR
22 mm,



SIMPLE STEM PLANE 5
JU-68-114
MOUTH MIRROR
24 mm,



SIMPLE STEM PLANE 4
JU-68-115
MOUTH MIRROR
22 mm,



SIMPLE STEM PLANE 4
JU-68-116
MOUTH MIRROR
22 mm,



SIMPLE STEM PLANE 4
JU-68-117
MOUTH MIRROR
22 mm,



SECO / CONE-SOCKET
JU-68-118
Mirror Handles
125 mm,5,

Mirror Handles



OCTOGONAL CONE SOCKET

JU-68-119

Mirror Handles
120 mm, 4,



MARTIN ERGO

JU-68-120

Mirror Handles
135 mm, 5?,



Mirror Handles

JU-68-121

Mirror Handles
135 mm, 5?,



Mirror Handles

JU-68-122

Mirror Handle
135 mm, 5?,



Mirror Handles

JU-68-123

Mirror Handles
135 mm, 5?,



Mirror Handles

JU-68-124

Mirror Handles
135 mm, 5?,

Mirror Handles


laryngeal Mirrors
JU-68-128
 MIRROR HANDLES
 7 mm 30°, ø 4,



laryngeal Mirrors
JU-68-129
 MIRROR HANDLES
 7 mm 30°, ø 4,



Cotton Wool Container
JU-68-130
 Cotton Wool Container
 ø 50 mm,



Napkin Holder
JU-68-131
 Napkin Holder
 ø 50 mm,



MIRROR HANDLE CONE SOCKET
JU-53-721-01
 MIRROR HANDLES
 14 cm /5,



MIRROR HANDLE CONE SOCKET
JU-53-721-02
 MIRROR HANDLES
 14 cm /5,

Mirror Handles



MIRROR HANDLE CONE SOCKET
JU-53-723-01
 MIRROR HANDLES
 13 cm /5,



MIRROR HANDLE
JU-53-728-02
 MIRROR HANDLES
 12 cm /4,



MARTIN-GRIP
JU-53-746-02
 12 cm /4?,



MIRROR HANDLE
JU-53-755-02
 MIRROR HANDLES
 14 cm /5 ø 11,



MIRROR HANDLE
JU-53-757-08
 MIRROR HANDLES
 14 cm /5,



MIRROR
JU-53-702-22

Mirror Handles



MIRROR
JU-53-703-22



RHODIUM MIRROR
JU-53-704-10
MIRROR



RHODIUM MIRROR
JU-53-704-18
MIRROR



CONE SOCKET
JU-53-706-22
MIRROR



CONE SOCKET
JU-53-707-22
MIRROR



MICRO-MOUTH MIRROR
JU-53-710-03

Fig. 1,


Mirror Handles



MICRO-MOUTH MIRROR

JU-53-710-05

Fig. 2,



RHODIUM MICRO MIRROR

JU-53-711-03

Fig.1 ? 3,0 mm,



HILGER

JU-53-715-00

18 cm /7,



HILGER

JU-53-715-01

CHEEK RETRACTORS

Size 1,

Endodontic Instruments



Fig. 1
Spreader

Root Spreader fig 1
JU-69-101
Endodontic Instruments
Fig. 1,



Fig. 2
Spreader

Root Spreader fig 2
JU-69-102
Endodontic Instruments
Fig- 2,



Fig. 3
Spreader

Root Spreader Fig 3
JU-69-103
Endodontic Instruments
Fig- 3,



Fig. 4
Spreader

Root Spreader Fig 4
JU-69-104
Endodontic Instruments
Fig- 4,



Fig. 5
Spreader

Root Spreader Fig 5
JU-69-105
Endodontic Instruments
Fig- 5,



Fig. 6
Spreader

Root Spreader Fig 6
JU-69-106
Endodontic Instruments
Fig- 6,

Endodontic Instruments



Fig. 7
Spreader

Root Spreader Fig 7
JU-69-107

Endodontic Instruments
Fig- 7,



Fig. 1
Luks

Luks Fig 1
JU-69-108

Endodontic Instruments
Fig- 1,



Fig. 2

LUKS Fig 2
JU-69-109

Endodontic Instruments
Fig- 2,



Fig. 3

LUKS Fig 3
JU-69-110

Endodontic Instruments
Fig- 3,



Fig. 4

LUKS Fig 4
JU-69-111

Endodontic Instruments
Fig- 4,



Fig. 1

Vertical Condenser Fig 1
JU-69-112

Endodontic Instruments
Fig- 1,

Endodontic Instruments



Fig. 2

Vertical Condenser Fig 2
JU-69-113
 Endodontic Instruments
 Fig- 2,



Fig. 3

Vertical Condenser Fig 3
JU-69-114
 Endodontic Instruments
 Fig- 3,



Fig. 6/1

Root Canal Explorer Fig 6/1
JU-69-116
 Endodontic Instruments
 Fig- 6/1,



Fig. 6/17

Root Canal Explorer Fig 6/17
JU-69-117
 Endodontic Instruments
 Fig- 6/17,



Fig. 6/23

Root Canal Explorer Fig 6/23
JU-69-118
 Endodontic Instruments
 Fig- 6/23,



Fig. 31L

Root Canal Excavators Fig 31L
JU-69-119
 Endodontic Instruments
 Fig- 31L,

Endodontic Instruments



Fig. 32L

Root Canal Excavators Fig 32L

JU-69-120

Endodontic Instruments
Fig- 32L,



Fig. 33L

Root Canal Excavators Fig 33L

JU-69-121

Endodontic Instruments
Fig- 33L,



Fig. 6/1

Root Canal Explorers Fig 6/1

JU-69-122

Endodontic Instruments
Fig- 6/1,



Fig. 6/17

Root Canal Explorers Fig 6/17

JU-69-123

Endodontic Instruments
Fig- 6/17,



Fig. 6/23

Root Canal Explorers Fig 6/23

JU-69-124

Endodontic Instruments
Fig- 6/23,



Fig. 31L

Root Canal Excavators Fig 31L

JU-69-125

Endodontic Instruments
Fig- 31L,

Endodontic Instruments



Fig. 32L

Root Canal Excavators Fig 32L

JU-69-126

Endodontic Instruments
Fig- 32L,



Fig. 33L

Root Canal Excavators Fig 33L

JU-69-127

Endodontic Instruments
Fig- 33L,

Separating Spatulas



**Heidemann
JU-70-101**
Separating Spatulas
Fig- 1,



**Heidemann
JU-70-102**
Separating Spatulas
Fig- 2,



**Heidemann
JU-70-104**
Separating Spatulas
Fig- 4,



**Heidemann
JU-70-105**
Separating Spatulas
Fig- 5,



**Heidemann
JU-70-106**
Separating Spatulas
Fig- 6,



**Heidemann
JU-70-107**
Separating Spatulas
Fig- 7,

Separating Spatulas



Heidemann
JU-70-108
Separating Spatulas
Fig- 8,



Heidemann
JU-70-109
Separating Spatulas
Fig- 9,



Heidemann
JU-70-110
Separating Spatulas
Fig- 1,



Heidemann
JU-70-111
Separating Spatulas
Fig- 2,



Heidemann
JU-70-112
Separating Spatulas
Fig- 3,

Excavators



Fig. 0

Excavators Fig 0
JU-71-101
Excavators
Fig- 0,



Fig. 1

Excavators Fig 1
JU-71-102
Excavators
Fig- 1,



Fig. 2

Excavators Fig 2
JU-71-103
Excavators
Fig- 2,



Fig. 3

Excavators Fig 3
JU-71-104
Excavators
Fig- 3,



Fig. 4

Excavators Fig 4
JU-71-105
Excavators
Fig- 4,



Fig. 5

Excavators Fig 5
JU-71-106
Excavators
Fig- 5,

Excavators



Fig. 6

Excavators Fig 6

JU-71-107
Excavators
Fig- 6,



Fig. 0

Excavators Fig 0

JU-71-108
Excavators
Fig- 0,



Fig. 1

Excavators Fig 1

JU-71-109
Excavators
Fig- 1,



Fig. 2

Excavators Fig 2

JU-71-110
Excavators
Fig- 2,



Fig. 3

Excavators Fig 3

JU-71-111
Excavators
Fig- 3,



Fig. 4

Excavators Fig 4

JU-71-112
Excavators
Fig- 4,

Excavators



Fig. 5

Excavators Fig 5

JU-71-113

Excavators
Fig- 5,



Fig. 6

Excavators Fig 6

JU-71-114

Excavators
Fig- 6,



Fig. 6

Excavators Fig 6

JU-71-115

Excavators
Fig- 6,



Fig. 14

Excavators Fig 14

JU-71-116

Excavators
Fig- 14,



Fig. 17

Excavators Fig 17

JU-71-117

Excavators
Fig- 17,



Fig. 17W

Excavators Fig 17W

JU-71-118

Excavators
Fig- 17W,

Excavators


Fig. 18

Excavators Fig 18

JU-71-119
Excavators
Fig- 18,



Fig. 18W
Oval spoon

Oval spoon Fig 18W

JU-71-120
Excavators
Fig- 18W,



Fig. 19

Excavators Fig 19

JU-71-121
Excavators
Fig- 19,



Fig. 19W
Oval spoon

Oval spoon Fig 19W

JU-71-122
Excavators
Fig- 19W,



Fig. 6

Excavators Fig 6

JU-71-123
Excavators
Fig- 6,



Fig. 14

Excavators Fig 14

JU-71-124
Excavators
Fig- 14,

Excavators



Fig. 17

Excavators Fig 17
JU-71-125
 Excavators
 Fig- 17,



Fig. 17W
 Oval spoon

Oval spoon Fig 17W
JU-71-126
 Excavators
 Fig- 17W,



Fig. 18W
 Oval spoon

Oval spoon 18W
JU-71-127
 Excavators
 Fig- 18W,



Fig. 19

Excavators Fig 19
JU-71-129
 Excavators
 Fig- 19,



Fig. 19W
 Oval spoon

Oval spoon Fig 19W
JU-71-130
 Excavators
 Fig- 19W,



Fig. 38/39

Excavators Fig 38/39
JU-71-131
 Excavators
 Fig- 38/39,

Excavators



Fig. 63/64

Excavators Fig 63/64

JU-71-133
Excavators
Fig- 63/64,



Fig. 65/66

Excavators Fig 65/66

JU-71-134
Excavators
Fig- 65/66,



Fig. 31L

Excavators Fig 31L

JU-71-135
Excavators
Fig- 31L,



Fig. 32L

Excavators Fig 32L

JU-71-136
Excavators
Fig- 32L,



Cleoid / Discoid

JU-71-137
No. 3/6 Excavator, double end



Cleoid / Discoid

JU-71-138
No. 89/92 Carver - Discoid \"/> Cleoid,
octagonal handle

Scalers


Fig. U15
Towner

Towner Fig U15
JU-72-104
Scalers
Fig- U15,



Fig. U15/33
Towner Jaquette

Towner Jaquette Fig- U15/33
JU-72-105
Scalers
Fig- U15/33,



Fig. U15/30
Towner Jaquette

Towner Jaquette Fig- 15/30
JU-72-106
Scalers
Fig- 15/30,



Fig H5/33

Fig- H5/33
JU-72-107
Scalers
Fig- H5/33,



Fig H6/H7

Fig- H6/H7
JU-72-108
Scalers
Fig- H6/H7,



Fig- 1/2

Fig- 1/2
JU-72-109
Scalers
Fig- 1/2,

Scalers



Fig. U15/33
Towner Jaquette

Towner Jaquette Fig- U15/33
JU-72-110
Scalers
Fig- U15/33,



Fig. J15/30
Towner Jaquette

Towner Jaquette Fig- U15/30
JU-72-111
Scalers
Fig- U15/30,




Fig. 1

Fig- 1
JU-72-112
Scalers
Fig- 1,



Fig. 33

JU-72-117
Scalers
Fig- 33, Fig- 31, Fig- 32, Fig- 34, Fig- 35,



Fig. H5
Hygienist

Hygienist Fig- H5
JU-72-122
Scalers
Fig- H5,

Wax & Modelling Carvers



P.K. THOMAS CARVERS FIG. 1

JU-73-101

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.1



P.K. THOMAS CARVERS FIG. 2

JU-73-102

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.2



P.K. THOMAS CARVERS FIG. 3

JU-73-103

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.3



P.K. THOMAS CARVERS FIG. 4

JU-73-104

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.4



P.K. THOMAS CARVERS FIG. 5

JU-73-105

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.5



P.K.THOMAS INSTRUMENTS

JU-73-106

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.5 SILVER

Wax & Modelling Carvers



**P.K.THOMAS INSTRUMENTS
JU-73-107**

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.4 RED



**P.K.THOMAS INSTRUMENTS
JU-73-108**

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.3 BLUE



**P.K.THOMAS INSTRUMENTS
JU-73-109**

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.2 GOLD



**P.K.THOMAS INSTRUMENTS
JU-73-110**

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.1 GREEN



JU-73-111

FAHNENSTOCK Wax Carver \\\\\" Knife,
5\\\\\\" , small



JU-73-112

FAHNENSTOCK Wax Carver \\\\\" Knife, 6-
3/4\\\\\\" , large

Impression Trays



Impression Trays
JU-74-101
Impression Trays
66x52 mm,



Impression Trays
JU-74-102
Impression Trays
62x50 mm,



Impression Trays
JU-74-103
Impression Trays
64x50 mm,



Impression Trays
JU-74-104
Impression Trays
65x46 mm,



Impression Trays
JU-74-105
Impression Trays
69x54 mm,



Impression Trays
JU-74-106
Impression Trays
65x46 mm,

Impression Trays



Impression Trays
JU-74-107
Impression Trays
70x48 mm,




Impression Trays
JU-74-108
Impression Trays
70x55 mm,



Impression Trays
JU-74-109
Impression Trays
59x42 mm,



Impression Trays
JU-74-110
Impression Trays
55x37 mm,



Impression Trays
JU-74-111
Impression Trays
61x42 mm,



Impression Trays
JU-74-112
Impression Trays
44 mm,

Impression Trays


Impression Trays
JU-74-113
 Impression Trays
 44 mm,



Impression Trays
JU-74-114
 Impression Trays
 44 mm,



Impression Trays
JU-74-115
 Impression Trays
 66x52 mm,



Impression Trays
JU-74-116
 Impression Trays
 65x46 mm,



Impression Trays
JU-74-117
 Impression Trays
 59x42 mm,



Impression Trays
JU-74-118
 Impression Trays
 69x54 mm,

Impression Trays



Impression Trays
JU-74-119
Impression Trays
70x48 mm,



Impression Trays
JU-74-120
Impression Trays
70x55 mm,



Impression Trays
JU-74-121
Impression Trays
55x37 mm,



Impression Trays
JU-74-122
Impression Trays
61x42 mm,



Impression Trays
JU-74-123
Impression Trays
61x42 mm,



Impression Trays
JU-74-124
Impression Trays
61x42 mm,

Impression Trays



Impression Trays
JU-74-125
Impression Trays

Dental Syringes



JU-75-101
Dental Syringes
55 x 37 mm,



JU-75-102
Dental Syringes
55 x 37 mm,




Self Aspirating 1.8cc
JU-75-103
Dental Syringes
55 x 37 mm,



Cylinder Cartridge Syringe
JU-75-104
Dental Syringes
55 x 37 mm,



JU-75-105
Dental Syringes
55 x 37 mm,



JU-75-106
Dental Syringes
55 x 37 mm,

Dental Syringes


JU-75-107
Dental Syringes
55 x 37 mm,



C-W 1.8cc
JU-75-108
Dental Syringes
55 x 37 mm,



JU-75-109
Dental Syringes
55 x 37 mm,



JU-75-110
Dental Syringes
55 x 37 mm,




JU-75-111
Dental Syringes
55 x 37 mm,



JU-75-112
Dental Syringes
55 x 37 mm,

Dental Syringes


JU-75-113
Dental Syringes
55 x 37 mm,



JU-75-118
Dental Syringes
Fig .1, Fig .2, Fig .3,



Dental Syringes
JU-75-119
Aspira-Plus
1.8 ml,



Dental Syringes
JU-75-120
Dental Syringes
1.8 ml,



Dental Syringes
JU-75-122
Dental Syringes
1,8 ml,



Dental Syringes
JU-75-124
Dental Syringes
1,8 ml,

Dental Syringes



Dental Syringes
JU-75-126
 Dental Syringes
 1,8 ml,



Dental Syringes
JU-75-127
 Dental Syringes
 1,8 ml,



Dental Syringe
JU-75-128
 Aspira-Plus
 1.8 ml,



Dental Syringe
JU-75-129



Dental Syringe
JU-75-130



Dental Syringe
JU-75-131

Dental Syringes



Dental Syringe
JU-75-132
1 ml,


Cavity Preparation

Wedeistaet
JU-76-101
Cavity Preparation

Fig- 1, Fig- 2, Fig- 3, Fig- 4, Fig- 1/2, Fig- 3/4, Fig- 5/6,



Wedeistaet
JU-76-112
Cavity Preparation

Fig- 1, Fig- 2, Fig- 3, Fig- 4, Fig- 1/2, Fig- 3/4, Fig- 5/6,

Explorers & Probes



11 Briault D/E Explorer
JU-77-101

FIG. # 11 BRIAULT D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



11/12 Odu D/E Explorer
JU-77-102

FIG. # 11/12 OLD DOMINION UNIVERSITY D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



11/12 Odu D/E After Five Explorer
JU-77-103

FIG. # 11/12 ODU D/E AFTER FIVE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



2A Pigtail D/E Explorer
JU-77-104

FIG. # 2A PIGTAIL D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



2H D/E Explorer
JU-77-105

FIG. # 2H DE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



2R/2L D/E Explorer
JU-77-106

FIG. # 2R/2L D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE

Explorers & Probes



3 Pigtail D/E Explorer
JU-77-107

FIG. # 3 PIGTAIL D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



3CH Cowhorn D/E Explorer
JU-77-108

FIG. # 3CH COWHORN D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



3ES D/E Explorer
JU-77-109

FIG. # 3ES DE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



5 D/E Explorer
JU-77-110

FIG. # 5 D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



6/3 D/E Explorer
JU-77-111

FIG. # 6/3 DE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



17 S/E Explorer
JU-77-112

FIG. # 17 SE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE

Explorers & Probes



23 Shepherd's Hook Explorer
JU-77-113

FIG. # 23 SHEPHERD'S HOOK S/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



3 S/E Explorer
JU-77-114

FIG. # 3 SE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



3A S/E Explorer
JU-77-115

FIG. # 3A S/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



54 S/E Explorer
JU-77-116

FIG. # 54 S/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE

Crown Removers



JU-78-101

A tips for MORREL Crown Remover



JU-78-102

B tips for MORREL Crown Remover



JU-78-103

C tips for MORREL Crown Remover



JU-78-104

MORREL Crown Remover, with 3 points



JU-78-105

WYMAN Carbide Crown Gripper, 5-1/2", without ratchet, with TC jaws

X-Rays Holders



X-Ray Film Hangers 6 Clips

JU-79-101

X-RAY FILM HANGERS 6 CLIPS MADE OF
STAINLESS STEEL



X-Ray Film Hangers 10 Clips

JU-79-102

X-RAY FILM HANGERS 10 CLIPS MADE
OF STAINLESS STEEL

Matrix Retainers


**Matrix Retainer Toffilemire
JU-80-101**

MATRIX BAND RETAINER TOFFILEMIRE.
AVAILABLE IN STRAIGHT & CONTRA.
STAINLESS STEEL



**Toffilemire Matrix Retainer Junior / Child
JU-80-102**

MATRIX BAND RETAINER TOFFILEMIRE
JUNIOR / CHILD. AVAILABLE IN
STRAIGHT & CONTRA. STAINLESS
STEEL



**Toffilemire Matrix Retainer Senior / Adult
JU-80-103**

MATRIX BAND RETAINER TOFFILEMIRE
SENIOR / CHILD. AVAILABLE IN
STRAIGHT & CONTRA. STAINLESS
STEEL



**Toffilemire Matrix Retainer Universal
JU-80-104**

MATRIX BAND RETAINER TOFFILEMIRE
UNIVERSAL. AVAILABLE IN STRAIGHT &
CONTRA. STAINLESS STEEL



**Matrix Band Retainer Ivory 8
JU-80-105**

MATRIX BAND RETAINER IVORY 8,
AVAILABLE IN 5MM, 6MM, 7MM.
STAINLESS STEEL



**Matrix Band Retainer Ivory Plain
JU-80-106**

MATRIX RETAINER IVORY PLAIN /
WITHOUT LATERAL SHIELD (JAW).
STAINLESS STEEL

Matrix Retainers



Matrix Retainer Ivory With Shield
JU-80-107

MATRIX RETAINER IVORY WITH LATERAL SHIELD (JAW). STAINLESS STEEL



Matrix Retainer Cheek Protector
JU-80-108

MATRIX BAND RETAINER CHEEK PROTECTOR. STAINLESS STEEL



Matrix Band Retainer Nystrom Left
JU-80-109

MATRIX BAND RETAINER NYSTROM LEFT. STAINLESS STEEL



Matrix Band Retainer Nystrom Right
JU-80-110

MATRIX BAND RETAINER NYSTROM RIGHT. STAINLESS STEEL



Matrix Band Retainer Siqveland (Ash Narrow)
JU-80-111

MATRIX BAND RETAINER SIQVELAND (ASH NARROW). STAINLESS STEEL



Matrix Band Retainer Siqveland (Ash Wide)
JU-80-112

MATRIX BAND RETAINER SIQVELAND (ASH WIDE). STAINLESS STEEL

Matrix Retainers



**Matrix Band Retainer Ivory Didner
(Adjustable)
JU-80-113**
Matrix Band Retainer Ivory Didner
(Adjustable)



**Matrix Retainer Ivory Separator
JU-80-114**
MATRIX BAND RETAINER IVORY
SEPARATOR. STAINLESS STEEL




**Matrix Retainer Elliott Separator
JU-80-115**
MATRIX RETAINER ELLIOTE
SEPARATOR. STAINLESS STEEL



**Matrix Retainer Ivory 8 N
JU-80-116**
MATRIX BAND RETAINER IVORY 8N.
MADE OF STAINLESS STEEL



**Matrix Retainer Didner 8 P
JU-80-117**
MATRIX BAND RETAINER IVORY 8P.
MADE OF STAINLESS STEEL

Articulators


PLASTERLESS ARTICULATOR GALETTI (WHITE)
JU-81-101

HIGH TECH PLASTERLESS ARTICULATOR GALETTI WITH ADJUSTABLE SCREWS AT REARS & FRONT. SOLIDLY CONSTRUCTED OF METAL ALLOY CASTING. ALL ANATOMICAL MOVEMENTS POSSIBLE



Vibratory Polished Articulator Galetti
JU-81-102

GALETTI ARTICULATOR IS SOLIDLY CONSTRUCTED OF METAL ALLOY CASTING. EASY TO USE, ALLOWING QUICK FIXATION OF MODEL. IMMEDIATE ADJUSTMENT OF THE OCCLUSAL & CORRECTS IMULATION OF CENTRAL RELATION. ALL ANATOMICAL MOVEMENTS POSSIBLE. VIBRATORY POLISHED



ANATOMICAL ARTICULATOR (ASA)
JU-81-103

ANATOMICAL ARTICULATOR WITH BEARING ROD & INCISIVE PIN, 4 MOUNTING PLATES & SUPPORTING BACK-PIN FOR 45° VIEW, OPENING 125°, INTERCONDYLAR DISTANCE 110 MM - SOLIDLY CONSTRUCTED OF METAL ALLOY CASTING



UNIVERSAL ARTICULATOR
JU-81-104

ANATOMICAL ARTICULATOR, MADE OF ALUMINUM PRESSOFUSION OF ONLY 590/GR OF WEIGHT. REALIZED FOR EASY VISIBILITY TO THE LINGUAL SURFACE & FOR A PERFECT STABILITY ALSO IN RESERVED POSITION. THE RETURN TO THE HINGS WITH RECALL SPRING IS AUTOMATIC. OPENING 125°, WITH 4 MOUNTING PLATES (AVAILABLE



Gnatus Pro JP30 Articulator
JU-81-105

The Gnatus Pro JP30 Articulator offers the best performance with an excellent degree of accuracy. Its functionality makes easier both the procedures of taking and transferring the patient's records and the process of mounting the models. Records can be taken and transferred accurately and fast. Light and resistant, makes operation easier. Its symmetrical and centered opening and



ANATOMICAL ARTICULATOR
JU-81-106

ANATOMICAL ARTICULATOR WITH INCISIVE PIN, 4 MOUNTING PLATES AND SUPPORTING BACK-PIN FOR 45° VIEW, OPENING 125°, INTERCONDYLAR DISTANCE 110 MM. - AVAILABLE IN ASSORTED COLOR

Articulators




KEYSTONE BLACK ARTICULATOR WIRED
JU-81-107

KEYSTONE DELUXE UNIVERSAL TYPE-A ARTICULATOR. INTERCHANGEABLE, INCISAL PIN INCLUDED, LATERAL EXCURSIONS.



KEYSTONE DELUXE MAGNET ARTICULATOR
JU-81-108

KEYSTONE DELUXE MAGNETIC DENTURE ARTICULATOR WITH PLATES, INTERCHANGEABLE, INCISAL PIN INCLUDED. MAGNETIC SPLIT CASTING ENSURES EASY & EXACT SETTING, LATERAL EXCURSIONS



PROFESSIONAL OCCLUSOR ARTICULATOR
JU-81-109

PROFESSIONAL OCCLUSOR, VIBRATORY POLISHED ARTICULATOR - ALUMINUM, BRASS, STAINLESS STEEL




Magicart-1 Articulator
JU-81-110

LIGHTWEIGHT, HIGH-QUALITY MAGIC ART ARTICULATOR SECURELY HOLDS QUADRANT MODELS, CONSERVE MATERIALS & TIME! QUICKLY & EASILY MOUNT PRE-CAST DENTAL MODELS WITH THE HIGH QUALITY. WITH MAGICART, YOU CAN SECURELY MOUNT & REMOVE ANY SIZE QUADRANT MODEL CASTS WITH AN EXCLUSIVE SIDE GRIP. SERRATED C-CLAMPS LOCK MODELS TIGHTLY WITH A QUICK ADJUSTMENT &



FULL DENTURE CHROME PLATED ARTICULATOR (3 SLOT ROUND)
JU-81-111



MODEL 'M' HI TOP STEPHAN ARTICULATOR
JU-81-112

FULL DENTURE CHROMIUM PLATED ARTICULATOR. PROVIDES ADDITIONAL VERTICAL DIMENSION FOR ACCEPTANCE OF HIGHER MODELS. 49

Articulators


Chromium Plated Plain Round
JU-81-113
 Nickle Chromium Plated Plain Round
 Articulator



UNIVERSAL ARTICULATOR (YELLOW)
JU-81-114
 UNIVERSAL ARTICULATOR OPENING 125°
 WITH 4 MOUNTING PLATES - CONDYLAR
 INCLINATION FIXED AT 30°,
 INTERCONDYLAR DISTANCE AT 110MM,
 BENNET ANGLE FIXED AT 7°, DISTANCE
 BETWEEN MOUNTING PLATES 85 MM.
 SOLIDLY CONSTRUCTED OF METAL
 ALLOY CASTING. Available In Assorted
 Colors



FOSTER DENTURE ARTICULATOR
JU-81-115
 THIS ARTICULATOR IS EQUIPPED WITH
 AN ADJUSTABLE ANTERIOR PIN STOP &
 POSTERIOR SCREW STOP, AND
 CONSTRUCTED OF A SLIGHTLY HEAVIER
 GAUGE BRASS GIVING THE
 ARTICULATOR A MORE SOLID FEEL.



FORK CHROME PLATED ARTICULATOR
JU-81-116
 NICKLE CHROMIUM PLATED
 ARTICULATOR FORK PATTERN



HINGE ARTICULATOR FULL ARCH
JU-81-117
 AN ECONOMICAL & SIMPLE CHROMIUM
 PLATED HINGE ARTICULATOR FOR
 MOUNTING FULL ARCH CASES. LOW



**QUADRANT / CROWNED BRIDGE
 ARTICULATOR**
JU-81-118
 BRASS HINGE CROWNED BRIDGE
 ARTICULATOR FOR MOUNTING
 QUADRANT CASES. CROWN & BRIDGE
 ARTICULATOR IS POPULAR FOSTER # 57
 DESIGN, ELIMINATES THE NEED TO
 BREAK OFF MODELS

Articulators



UNIVERSAL ARTICULATOR (BLUE)
JU-81-119

UNIVERSAL ARTICULATOR OPENING 125° WITH 4 MOUNTING PLATES - CONDYLAR INCLINATION FIXED AT 30°, INTERCONDYLAR DISTANCE AT 110MM, BENNET ANGLE FIXED AT 7°, DISTANCE BETWEEN MOUNTING PLATES 85 MM. SOLIDLY CONSTRUCTED OF METAL ALLOY CASTING.



FOSTER / APEX # 2 ARTICULATOR
JU-81-120

FOSTER / APEX # 2 ARTICULATOR, THE WIDTH AT THE JOINT IS 1-3/4" WITH A VERTICAL PARALLEL HEIGHT OF 2-1/8". THE SMALL SIZE OF THIS ARTICULATOR ALLOWS THE VERSATILITY TO USE FOR FULL OR PARTIAL WORK. THE FINISH IS BRIGHT CHROME. HIGH BACK



MORTIMER-I OCULADER ARTICULATOR
JU-81-121

PLAINLINE BALL BEARING MORTIMER-I OCULADER ARTICULATOR. ALL FILLED WITH HIGH CARBON STEEL BALLS FOR SMOOTH SETTING & ARTICULATION OF DENTURE



MORTIMER-II OCULADER ARTICULATOR
JU-81-122

PLAINLINE BALL BEARING MORTIMER-II OCULADER ARTICULATOR. ALL FILLED WITH HIGH CARBON STEEL BALLS. SEMI-ADJUSTABLE ARTICULATOR ALLOWS THE ESTABLISHMENT OF MORE ACCURATE OR CUSTOMISED LATERAL & PROTRUSIVE AS WELL AS CENTRIC OCCLUSION



OCCLUSOR, WIRE ARTICULATOR
JU-81-123

OCCLUSOR, WIRE ARTICULATOR - STAINLESS STEEL & ALUMINUM



ASA Articulator 5002
JU-81-124

ASA ARTICULATOR WITH LATERAL MOMENT. MADE OF ALUMINUM DIE CASTING WITH 2 MOUNTING PLATES. ALSO AVAILABLE IN VIBRATORY POLISHED PROCESS



Model \"A\" Crescent Stephan Articulator

JU-81-125

CRESCENT STEPHAN MODEL \"A\" ARTICULATOR IS USED FOR PARTIALS, CENTER SPRING WITH ADJUSTABLE TENSION, SMOOTH BRASS CASTINGS WITH NICKEL PLATING, CENTER BAR IN BOTH UPPER & LOWER BOW FOR MAXIMUM PLASTER RETENTION & STRENGTH. WIDTH AT THE JOINT OF 29, & VERTICAL SPAN WHEN BOWS ARE PARALLEL OF 2 1/8\". REPRODUCE ALL ANATOMICAL MOVEMENTS. UPPER AND LOWER SECTIONS MAY BE SEPARATED OR ASSEMBLED QUICKLY. OPEN CONSTRUCTION PROVIDES FREE ACCESS TO LINGUAL SURFACES. ANATOMICAL IN EVERY RESPECT



Bite Fork / Occlusal Alignment Plane

JU-81-126

BITE FORK IS A RELIABLE INSTRUMENT FOR CHECKING THE CORRECT POSITION OF THE MASTICATORY PLANE, ACHIEVING PARALLELISM OF THE OCCLUSAL PLANE WITH THE CAMPER'S LINE & THE BIPUPILLARY LINE IS FUNDAMENTAL TO THE CONSTRUCTION OF DENTURES

Miscellaneous Laboratory Instruments



JU-82-101

Bandage and Utility Scissors, 7-1/2" (19.1 cm), needle destroyer, serrated blade, autoclavable



JU-82-102

Plaster and Compound Knife



JU-82-103

Soldering Tweezers, curved



JU-82-104

SWISS Jeweler style Forceps, non-magnetic stainless steel, style 3, narrow, fine, 4-3/4" (12.1 cm)



JU-82-105

SWISS Jeweler style Forceps, non-magnetic stainless steel, style 5, Super fine, 4-3/8" (11.2 cm)



JU-82-106

SWISS Jeweler style Forceps, non-magnetic stainless steel, style 7, fine, curved, 4-1/2" (11.4 cm)

Miscellaneous Laboratory Instruments



Calipers
JU-82-107
Calipers



Calipers
JU-82-108
Calipers



Steel Ruler
JU-82-109
Steel Ruler



Forceps for articulation paper
JU-82-110
Forceps for articulation paper



CONTACT US

Jimy Brothor Ins

Address:

339A North High Street MT Vernon
10550, USA

T: 347 469 5700

F: 347 469 5700

PRODUCTS CATALOGUE